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# PC15

## Phosphate Investment for Pressable Ceramics

Formulated for use with both staining and layering techniques

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
## Recommended Special Liquid: FastFire Liquid

## 1. Mixing:

- Add measured liquid and weighed powder to bowl. Hand incorporate 15–20 seconds, mechanical mix under vacuum, slow speed, for 45 seconds (350-600 rpm).
- Vibrate mix lightly under vacuum (5 seconds) then slowly release vacuum.

## 2. Investing:

- Slowly pour investment into ring using low vibration. Follow glass manufacturer's recommendations for filling and preparing invested ring.  
**Note: To avoid trapping air bubbles, use a sable brush or small blunt instrument to paint and fill wax patterns prior to filling ring. The use of a wetting agent/debubblizer is not recommended.**

## 3. Benchset:

- Benchset a **minimum** of 15 minutes.  
**If invested ring(s) are allowed to benchset over a 12-hour period, re-wet mold by soaking in water 3–5 minutes.**
- Carefully remove leveling ring, sprue former, and mold from ring. Using a lab knife, level any excess nub of investment from top of mold (as a result of leveling ring) and carefully smooth seam line formed in investment mold by ring.

## 4. Burnout:

## Rapid Technique

- Place cylinder(s) in preheated oven at 850°C (1562°F).
- For systems requiring higher burnout temperature, place mold into preheated oven at 850°C (1562°F) — hold for 5 minutes and raise temperature to desired/ recommended temperature at 8°C per minute.
- Heat soak for 45 minutes for 100 g mold or 60 minutes for 200 g mold (add 10 minutes for each additional mold).

**Note: If 100% liquid concentration is required for expansion purposes, a lower preheat temperature is recommended. Place molds in preheated oven at 650°C (1200°F) to 850°C (1562°F) and heat soak at 14–20°C per minute.**

## Standard Technique:

- Place in cold oven.
- Heat to 850°C (1562°F) at 14–20°C (25–40°F) per minute.
- Heat soak for 45 minutes for 100g mold or 60 minutes for 200g mold (add 10 minutes for each additional mold).

## 5. Pressing:

- Process glass according to manufacturer's recommendations.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

## Physical Properties:

Liquid/Powder Ratio	27 mL/100 g
Initial Setting Time	7–9 minutes
Setting Expansion	1.0%
Thermal Expansion	0.7%
Compressive Strength, Wet (45 minutes)	500 psi (3.4 MPa)

\*Measured at 100% liquid concentration.

## Suggested FastFire Liquid Concentrations\*

Restoration	Concentration	100 g Investment/27 mL Special Liquid: Distilled Water	200 g Investment/54 mL Special Liquid: Distilled Water
Crowns/Veneers	80%	22 mL liquid/5 mL water	44 mL liquid/10 mL water
Inlays — MO, DO	60%	16 mL liquid/11 mL water	32 mL liquid/22 mL water
Inlays — MOD	75%	20 mL liquid/7 mL water	40 mL liquid/14 mL water
3/4 Crowns, Onlays			

\* Recommended concentrations are approximate and can be adjusted as required. To increase expansion, use more liquid and less water. To decrease expansion, use less liquid and more water. Always maintain a total of 27 mL liquid/water volume per 100 grams. Do not dilute below 50%.

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