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FastFire 15

Universal Phosphate Investment for All Crown and Bridge Alloys

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
ENGLISH

1. Pattern Preparation:

- Attach sprue to the wax pattern using pliable sticky wax.
- Mount patterns on crucible former base.
- Spray with Whip Mix SMOOTHX Wax Pattern Cleaner. Gently blow excess cleaner from pattern.

2. Preparation:

- If metal ring is used, line casting ring with a ceramic or paper liner positioned 6 mm short of each end.
- Soak paper lined ring in water for a minimum of one minute, then shake excess from ring before investing (DO NOT WET THE CERAMIC LINER).
- Place ring onto former base.

3. Mixing:

Note: For best results store liquid and powder at temperatures of 20–25°C (68–78°F).

- Prepare liquid at suggested concentration following the chart below (distilled water recommended for dilution).
- Rinse bowl out with water and shake out excess. Use separate mixing bowls for phosphate and gypsum investments.
- Add measured liquid to mixing bowl. Incorporate powder by hand spatulation 10–15 seconds.
- Mechanical mix under vacuum, slow speed (350–600 RPM) for **60 seconds**. Higher RPM mixers may require decreased mix time.

Note: Different paddle designs and mixer speed may require varying mixing time.

4. Benchset:

- Benchset 15 minutes minimum.
Note: When using 60 grams or less, benchset for 30 minutes.
- Trim glaze off top of mold and rinse under tap water.
- For best results, place in preheated oven within 20 minutes of investing.

Note: Molds allowed to set more than 12 hours should be re-wet by soaking in water for 1–3 minutes prior to burnout.

5. Burnout:

Note: If using ringless casting system, remove mold at 15 minutes, then allow to benchset an additional 5 minutes for 100g mold and 10 minutes for 200g mold after removal from the ring. Use rapid burnout method for benchset up to 1 hour. For longer than 1 hour benchset, conventional burnout is recommended.

Rapid Technique (Preheated Oven):

- Place molds in preheated oven at alloy manufacturer's recommended temperature, up to 900°C (1650°F). For higher temperatures, place molds in oven at 900°C (1650°F) then heat to final temperature at 14–20°C (25–35°F) per minute.
- Heat soak at final temperature for 30 minutes, add 10 minutes per additional mold.

Standard Technique (Cold Oven):

- Place molds in oven at room temperature.
- Heat to desired temperature at 14–20°C (25–35°F) per minute.
- Heat soak at final temperature for 30 minutes, add 10 minutes per additional mold.

6. Casting:

- Upon removal from oven, immediately cast according to alloy manufacturer's instructions.
- Use separate mixing bowls for phosphate and gypsum investments.
- To correct tight fits: Increase liquid concentration. If using metal rings you may also use a double liner.
- To correct loose fits: Decrease liquid concentration.
- When using bulk packaging, reseal container to avoid physical property change.
- For large molds containing complex restorations or plastic sprues, runner bars or copings, the standard technique described above or a two-stage burnout technique is recommended.

WARNING: Investments contain free silica — DO NOT BREATHE DUST. May cause delayed lung injury (silicosis/lung cancer).

Physical Properties:

Liquid/Powder Ratio	27 mL/100 gram
Working Time	7–9 minutes
Setting Expansion	1.0%
Thermal Expansion	1.0%
Compressive Strength	500 psi (3.4 MPa)

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Suggested Liquid Concentrations*

Alloy	Liquid Concentration	16 mL/60 g	24 mL/90 g	27 mL/100 g
Base (Ni-Cr, Ni-Cr-Be)	75%	12 mL liquid 4 mL water	18 mL liquid 6 mL water	21 mL liquid 6 mL water
Noble (Ag-Pd, Au < 40%)	65%	10 mL liquid 6 mL water	15 mL liquid 9 mL water	17 mL liquid 10 mL water
High Noble (Pd w/ 2–8% Au, Au > 40%)	50%	8 mL liquid 8 mL water	12 mL liquid 12 mL water	14 mL liquid 13 mL water

* Recommended concentrations are approximate and can be adjusted to optimize fit. To increase expansion, use more liquid and less water. To decrease expansion, use less liquid and more water. Always maintain total liquid/water volume.

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