## UNIFAST™ LC

LIGHT CURED ACRYLIC RESIN

For use only by a dental professional in the recommended indications.

### INDICATIONS FOR USE

Direct technique using a hand held visible light curing device (470 nm wavelength).

Indirect technique using a table top curing unit (470 nm wavelength).

By pressure and/or brush on techniques.


- Temporary crowns, bridges, inlays and onlays.
- Fixing the position of metal cast crowns for soldering.
- Repair of dentures.
- Individual tooth trays, individual trays.
- Orthodontic plates.
- Check bite, splint, rebase, night guard, etc.

### CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

### DIRECTIONS FOR USE

- The standard powder/liquid ratio is 1.0 g/0.5 mL (first graduation on the powder measure/21 drops of liquid).
- Add powder to liquid in the mixing cup and mix quickly for 10-15 seconds.
- The following illustrates the lapsed time stages in polymerisation:

- Use the slurry or dough according to application standard techniques.
- When all necessary corrections or adjustments have been made, light cure the resin using an appropriate visible light curing device (470 nm wavelength).

Note:

- Quality of cure, and depth of cure, vary depending on the light source used, and on the translucency or opacity of the shade of resin.

|                       | Depth of cure (mm) |       |            |               |
|-----------------------|--------------------|-------|------------|---------------|
|                       | Translucent        | Ivory | Pink # 3/E | Live Pink # 8 |
| Hand held: 20 seconds | 6-9                | 3-4   | 2-3        | 3-4           |
| 60 seconds            | 13-16              | 5-6   | 3-4        | 5-7           |
| Table top: 1 minute   | 10                 | 3     | 2          | 5             |
| 10 minutes            | >25                | 16    | 16         | 24            |

The above are typical values and are given for guidance only.

- Light cure within 1 hour from starting the mix.
- UNIFAST LC can also be cured chemically in 5 minutes in water at 60°C (140°F).

### STORAGE

Recommended for optimal performance, store in a cool dark place away from direct sunlight (4-25°C) (39-2-77.0°F).

### PACKAGES

- Bottle of 50 g powder.
- Box of 30 g liquid (2 x 15 g or 14.3 mL).
- 1-2 package\*: 1 x 50 g powder, 2 x 15 g (14.3 mL) liquid, accessories.

\* Not available in Europe.

### SHADES

- Crown colour Ivory\*, A2, A3, B2, B3, C2
- Gingival colour\* Pink # 3, Live Pink # 8 (veined), Pink E
- Translucent Translucent

\* Not available in Europe.

### CAUTION

- The liquid is volatile and highly flammable. Close the bottle immediately after use and keep away from sources of ignition.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Do not use powder or liquid with components of other brands of acrylic resin.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people.
- If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarise yourself with the safety data sheets available at:

<http://www.gceurope.com>

or for The Americas

<http://www.gcamerica.com>

They can also be obtained from your supplier.

### CLEANING AND DISINFECTING

MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

For the Summary of Safety and Clinical Performance (SSCP) please see EUDAMED database (<https://ec.europa.eu/tools/eudamed>) or contact us at [Regulatory.gce@gc-dental](mailto:Regulatory.gce@gc-dental)

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

[https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en)

as well as to our internal vigilance system: [vigilance@gc-dental](mailto:vigilance@gc-dental)

In this way you will contribute to improve the safety of this product.

UK Responsible Person  
GC UNITED KINGDOM Ltd.  
Coopers Court Newport Pagnell, Buckinghamshire, MK16 8JS, United Kingdom  
Last revised: 10/2022

**RxOnly**

**IGC** 2797 **UK** 0086

MANUFACTURED BY

GC DENTAL PRODUCTS CORP.

2-095 Tomimachi-cho, Nakagyo, Aichi 464-0844, Japan

RESPONSIBLE MANUFACTURER IN CANADA:

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