


# EXA'lence™ 48mL User Guide



Bleed cartridge before first use.



Make sure to have no excess materials at the end of the cartridge.



Store with the used mixing tip on.

GC has harnessed the best qualities of two great materials to form the next generation of impression material: Vinyl PolyEther Silicone (VPES™).

Because of the unique chemistry, we provide this User Guide for some special handling tips.

EXA'lence clinical benefits include: intrinsically hydrophilic composition, superior flow in the sulcus, superior elasticity without deformation, high tear strength and pleasant mint taste.

## Mixing Tip Guide:

- For Extra Light Body and Light Body use Yellow Tips
- For Monophase, Heavy Body and Heavy Body Rigid use Aqua Tips



# EXA'lence™ 48mL User Guide

**Total working time:**  
**Regular set: 2 min.**  
**Fast set: 1 min.**

**Minimum time in mouth:**  
**Regular set: 3 min.**  
**Fast set: 1.5 min.**

**Note:**

EXA'lence has unique chemistry that allows it to polymerize faster in the oral cavity than at room temperature. Therefore to determine setting, examine material in the mouth since material on bench will yield a false result.

