Miracle Mix™ CAPSULE

RADIOPAQUE SILVER ALLOY-GLASS IONOMER CEMENT RESTORATIVE IN CAPSULES

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

1. Core build-up.
2. Restoration of primary teeth.
3. Cases where radiopacity is required.
4. As a base or liner.
5. Class I, limited Class II, temporary fillings, etc.

CONTRAINDICATIONS

1. Pulp capping.
2. In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

DIRECTIONS FOR USE

1. CAVITY PREPARATION

- a) Prepare tooth using standard techniques. Extensive mechanical retention is not necessary. For pulp capping, use calcium hydroxide.
- b) Wash and dry but DO NOT DESICCATE. GC CAVITY CONDITIONER (10 seconds) or GC DENTIN CONDITIONER (20 seconds) is recommended to remove the smear layer (Fig. 1).

2. CAPSULE ACTIVATION / MIXING

- a) Before activation shake the capsule or tap its side on a hard surface to loosen the powder (Fig. 2).
- b) To activate the capsule, push the plunger until it is flush with the main body and hold it down for 2 seconds (Fig. 3).
Nota:
 - i) Ensure the plunger is fully pressed to avoid the incorrect mixing ratio of powder and liquid.
 - ii) The capsule should be activated just before mixing and used immediately.
- c) Immediately remove the capsule and set it into a mixer (an amalgamator) and mix for 10 seconds at high speed (~4,000 RPM) (Fig. 4).

3. RESTORATIVE TECHNIQUE

- a) Immediately remove the mixed capsule from the mixer and load it into the GC Capsule Applier. Time of initial extrusion after mixing : 10 seconds maximum.
- b) Make two clicks to prime the capsule then syringe (Fig. 5). The working time is 1 minute from start of mixing at 23°C (73.4°F). Higher temperature will shorten working time.
- c) Remove surface moisture but DO NOT DESICCATE.
- d) Within 10 seconds maximum after mixing, start to extrude the mixture directly into the preparation (Fig. 6). Avoid air bubbles.
 - 1) To adjust the direction of the nozzle, hold the applicator with the capsule towards you and turn the capsule body.
 - 2) To remove the used capsule, push the applicator release button. Twist the capsule and pull upwards.
- e) Form the contour.
- f) Setting time is 5'00" from start of mixing.

4. FINISHING

- a) Use standard techniques (Fig. 7).

STORAGE

Recommended for optimal performance, store the Miracle Mix capsules in the original aluminum foil in a cool and dark place (4-25°C)(39.2-77.0°F).

PACKAGE

Package: 48 capsules.

Average contents per capsule: 0.52g silver alloy ionomer cement and 0.11g (0.09mL) liquid.
Minimum net volume of mixed cement per capsule: 0.09mL.
Single Shade Export package: 100 capsules
Option: GC Capsule Applier (1 piece).

CAUTION

1. In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water.
2. In case of contact with eyes, flush immediately with water and seek medical attention.
3. Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available for the Americas at: <http://www.gcamerica.com>
They can also be obtained from your supplier.

Last revised: 11/2019

GC

MANUFACTURED by
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-ku, Tokyo 174-8585, Japan

DISTRIBUTED by
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A.
www.gcamerica.com
TEL: +1-708-597-0900

PRINTED IN THE U.S.A.

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3737 West 127th Street, Alsip, IL 60803 U.S.A

GC SOUTH AMERICA
Rua Heliodora, 399, Santana - São Paulo, SP -
Brasil CEP: 02022-051 - TEL: +55-11-2925-0965
CNPJ: 08.279.999/0001-61
RESP. TEC. Mayara De Santos Ribeiro -
CROSP/ 105.982

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