
# **ceramill motion 2**

## **Gebrauchsanweisung Instruction Manual Instructions d'utilisation Istruzioni per l'uso Instrucciones de uso**

Deutsch	3 - 16
English	17 - 30
Français	31 - 44
Italiano	45 - 58
Español	59 - 72

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Tab. 3

- Translation of the original Operating Instructions -

## Table of Contents

<b>1</b>	<b>Product description</b>	<b>17</b>
1.1	Delivery scope	18
1.2	Required system components	18
<b>2</b>	<b>Components</b>	<b>18</b>
<b>3</b>	<b>Calibrating the blank holder</b>	<b>19</b>
3.1	Zero-point checking of the blank	19
3.2	Converting the Ceramill Motion 2 (5X) for calibration milling	19
3.2.1	Converting the machine for wet operation	19
3.2.2	Removing the blank holder from the machine	19
3.2.3	Inserting the Ceramill blank holder for rotational milling	20
3.3	Calibration process	21
3.3.1	Milling the measuring surfaces	22
3.3.2	Measuring the measuring surfaces	23
3.3.3	Performing a calibration milling	24
3.3.4	Mounting the Ceramill blank holder for rotational milling	26
3.3.5	Completing the calibration process	27
3.4	Changing a calibration sleeve	27
<b>4</b>	<b>Processing Ceramill Ti form blanks</b>	<b>28</b>
4.1	General workflow	28
4.2	Tool placement	28
4.3	Blank positioning	29
<b>5</b>	<b>Spare parts</b>	<b>30</b>

## 1 Product description

Calibratable blank holder for use in the Ceramill Motion 2 (5X).

The Ceramill blank holder for rotational milling (blank holder 4) consists of two components:

- \_ Tool holder including calibration sleeve
- \_ Blank holder

After a calibration milling process, both components are screwed together.

The blank holder 4 is used to affix prefabricated titanium abutment blanks (Ceramill Ti forms) and to take up appropriate tools (Ceramill Roto Ti).

Ceramill Ti form blanks are available for various implant systems (see catalogue and [www.amanngirrbach.com](http://www.amanngirrbach.com)).

In order to process Ceramill Ti form blanks, the machine must be converted and prepared for wet operation (see Ceramill Motion 2 Operating Instructions).

## 1.1 Delivery scope


Fig. 1

- 1 Tool holder including calibration sleeve
- 2 Additional calibration sleeve
- 3 Blank holder
- 4 Allen screws, M4
- 5 Allen key (size 3)
- 6 Allen key (size 2.5)

## 1.2 Required system components

Art. No.	Component
179156	Ceramill Mind/Match 2
179152	Ceramill M-Plant
179250	Ceramill Motion 2 (5X)
178630	Ceramill Coolstream
760615	Ceramill Roto TI 1.0
760616	Ceramill Roto TI 2.0
760301	Ceramill TEST 71 L
	Caliper (delivery scope of Ceramill Motion 2)

Tab. 1

In order to ensure the software functionality:

► Install the current update of the following

Ceramill software components:

- Ceramill Mind
- Ceramill Match 2
- Ceramill Motion

The current updates can be downloaded in the protected download area of the Amann Girrbach homepage:

<https://www.amanngirrbach.com/de/services/downloads/>

## 2 Components


Fig. 2 Components of the Ceramill blank holder for rotational milling

- 1 Connection mandrel
- 2 Tool holder
- 3 Calibration sleeve (mounted to tool holder)
- 4 Blank holder
- 5 Bolt for anti-rotation lock
- 6 Allen screws (M 4)

## 3 Calibrating the blank holder

Before using the Ceramill blank holder for rotational milling in the Ceramill Motion 2 (5X), the calibration sleeve must be calibrated. Afterwards, the blank holder is slid over the calibration sleeve and screwed to the tool holder.

To ensure correct calibration of the tool holder, the following steps must be carried out in the sequence described:

- ▷ Perform a zero-point check of the blank.
- ▷ convert the Ceramill Motion 2 (5X) for calibration milling.
- ▷ Perform the calibration milling.

### NOTE:

Insufficient milling results!

The calibration sleeve of the blank holder is individually adjusted to the respective milling machine.

- ▷ Do not use a calibrated Ceramill blank holder for rotational milling in other machines!

### 3.1 Zero-point checking of the blank

### NOTE:

Insufficient milling results!

Without checking the zero point, the calibration process will not be carried out correctly.

Before carrying out the calibration process:

- ▷ Check the zero point of the blank in the Ceramill Motion (see "Ceramill Calibration Process.pdf" in the protected download area of [www.amanngirrbach.com](http://www.amanngirrbach.com)).

### 3.2 Converting the Ceramill Motion 2 (5X) for calibration milling

#### 3.2.1 Converting the machine for wet operation

The milling machine must be equipped for wet operation (see Ceramill Motion 2 Operating Manual).

#### 3.2.2 Removing the blank holder from the machine

- ▷ Start the Ceramill Motion software and connect to the Ceramill Motion 2 (5X).
- ▷ Click the *Service and control functions*  button.

- ▷ Activate the *Deposit tool and change blankholder* function.

The current tool is placed down and the blank holder [2] moves to the change position.


Fig. 3 Removing the extraction hood

- 1 Extraction hood
- 2 Blank holder
- 3 Lever
- A Operation button

- ▷ Release lever [3] of extraction hood [1] and pull off the extraction hood toward the front.
- ▷ Loosen fastening screw [1] and pull off the blank holder 1 toward the right.


Fig. 4 Removing the blank holder 1

- 1 Fastening screw for blank holder
- 2 Blank holder 1

- ▷ Clean the locating bore and surfaces.

### 3.2.3 Inserting the Ceramill blank holder for rotational milling

- ▷ Check the connection mandrel of the blank holder and clean if required.
- ▷ Insert the Ceramill Roto TI tools into the Ceramill tool holder for rotational milling with the tool shank facing outward (see chapter 4.2 on page 28).
- ▷ Insert the tool holder together with the calibration sleeve via the connection mandrel to the stop into the locating bore.


Fig. 5

- ▷ Tighten fastening screw [1] with the torque wrench until the torque wrench clicks once.


Fig. 6 Fastening the tool holder 4

- 1 Fastening screw for tool holder
- 2 Tool holder 4

- ▷ Mount the transparent splash guard.
- ▷ Insert the swarf sieve.

► Select 4 (*TI-Forms*) in the Ceramill Motion Software and confirm with *OK*.


Fig. 7

The machine picks up the tool from location 1 and advances to the home position.

The machine is now prepared for the calibration-milling process.

## 3.3 Calibration process

The calibration process can only be performed once per calibration sleeve. Should errors occur during the process, we recommend replacing the calibration sleeve (see chapter 3.4 on page 27) and starting the process from the beginning.

The calibration process consists of 3 steps that must be performed one after the other.

1. Milling the measuring surfaces
2. Measuring the measuring surfaces (e.g. with the caliper in the delivery scope of the Ceramill Motion 2)
3. Performing the calibration milling

## 3.3.1 Milling the measuring surfaces

- ▷ Start the Ceramill Motion software and connect to the Ceramill Motion 2 (5X).
- ▷ Click on the *Open program file(s)*  button.
- ▷ Under the *Test Programs* tab, select and transfer the *Blank-Holder-4-Calibration\_1* program.


Fig. 8

- ▷ Make sure that the machine is set up according to the necessary data of the following window and confirm with *OK*.  
The milling data are transferred.
- ▷ After completion of the data transfer, reset the tool wear for all tools.

Fig. 9

- ▷ To start the milling process, click on the  button.  
Two plane surfaces are milled into the calibration sleeve.

The milling procedure is completed after approx. 20 minutes and the operation button [A] starts to flash. After the milling process, the tool holder remains in the removal position.

 There is **no** indication in the Ceramill Motion software that the process has been completed.

- ▷ Remove the tool holder from the machine.
  - Release the lever of the splash guard and pull off the splash guard toward the front.
  - Loosen the fastening screw of the tool holder and pull off the tool holder toward the right.
- ▷ Clean the locating bore and the surfaces of the machine.

## 3.3.2 Measuring the measuring surfaces


Fig. 10

Two plane and opposite surfaces have been milled in the rear third of the uncoated surface of the calibration sleeve. The distance between them must be measured.

### NOTE:

Faulty calibration!

► An accurate measurement is required for a successful calibration.

- Centrally position the external jaw to one of the two surfaces and measure the distance to the opposite surface with the moving internal jaw. The milled surfaces of the calibration sleeve must be fully embraced by both caliper jaws.
- Note down the determined measuring value two digits behind the decimal point, e.g., 31.93 mm.

- If required, clean the caliper jaws.
- Set the unit of measurement on the caliper to mm.
- Slide the caliper jaws together and set to 0.

## 3.3.3 Performing a calibration milling

- ▷ Check the connection mandrel of the Ceramill blank holder for rotation milling and clean if required.
- ▷ Insert tooling into tool holder.

The tool of location 1 is already in the collet.

- ▷ Insert the tool holder via the connection mandrel to the stop into the locating bore of the machine.


Fig. 11

- ▷ Tighten the fastening screw with the torque wrench until the torque wrench clicks once.
  - ▷ Mount the transparent splash guard to the extraction opening and lock it in place with the lever.
  - ▷ Briefly press the flashing operation button [A]. The machine advances to the home position.
  - ▷ In the Ceramill Motion software, click on the *Program successfully finished* button.
  - ▷ Click on the *Open program File(s)*  button.
  - ▷ Open the *Test Programs* tab.
- Several program files are listed.


Fig. 12

MM	FN	MM	FN	MM	FN	MM	FN
32.12	..._32.12	32.05	..._32.05	31.98	..._31.98	31.91	..._31.91
32.11	..._32.11	32.04	..._32.04	31.97	..._31.97	31.90	..._31.90
32.10	..._32.10	32.03	..._32.03	31.96	..._31.96	31.89	..._31.89
32.09	..._32.09	32.02	..._32.02	31.95	..._31.95	31.88	..._31.88
32.08	..._32.08	32.01	..._32.01	31.94	..._31.94		
32.07	..._32.07	32.00	..._32.00	31.93	..._31.93		
32.06	..._32.06	31.99	..._31.99	31.92	..._31.92		

Tab. 2 Assignment of measuring results - Program file

**MM** Measuring result in mm

**FN** End of file name (Blank-Holder-Calibration...)

► Select the program, whose last 4 file name digits correspond with the noted down measuring value in step 1. E.g.:

- Determined measuring value MM: 31.93 mm
- FN calibration file to be selected: Blank-Holder-Calibration\_31.93

► Transfer the respective program to the machine.

► Make sure that the machine is set up according to the necessary data of the program window and confirm with *OK*.

The milling data are transferred.

► After completion of the data transfer, reset the tool wear for all tools.

The milling procedure is completed after approx. 20 minutes and the operation button [A] starts to flash. After the milling process, the tool holder remains in the removal position.

There is **no** indication in the Ceramill Motion software that the process has been completed.

► Remove the tool holder from the machine.

- Release the lever of the splash guard and pull off the splash guard toward the front.
- Loosen the fastening screw of the tool holder and pull off the tool holder toward the right.

► Clean the locating bore and the surfaces of the machine.

The calibration milling process is completed.

Fig. 13

► To start the milling process, click on the button.

The necessary fitting surfaces for applying the blank holder are milled into the calibration sleeve.

## 3.3.4 Mounting the Ceramill blank holder for rotational milling

### NOTE:

When assembling, the anti-rotation bolt in the blank holder must be in the provided groove in the calibration sleeve.

▷ Avoid jamming/wedging!

▷ Carefully mount the blank holder onto the calibration sleeve of the tool holder and firmly screw together with the two provided screws using a size 3 Allen key.

After screwing the blank holder and the tool holder together, a small gap will remain, which is negligible.

The Ceramill blank holder for rotational milling can only be used for this machine.  
▷ Make an appropriate mark on the blank holder.


Fig. 14

## 3.3.5 Completing the calibration process

To complete the calibration process:

- ▷ Clamp the assembled blank holder 4 into the machine.
  - Check the connection mandrel of blank holder 4 and clean if required.
  - Insert blank holder 4 via the mandrel to the stop in the locating bore.
  - Tighten the fastening screw with the torque wrench.
- ▷ Mount the transparent splash guard to the extraction opening and lock it in place with the lever.
- ▷ Briefly press operation button [A].
 

The machine advances to the home position.
- ▷ In the Ceramill Motion software, click on the *Program successfully finished* button.
 

The blank holder 4 is now calibrated for this machine and can immediately be used for machining Ceramill Ti form blanks.

## 3.4 Changing a calibration sleeve

The Ceramill blank holder for rotational milling is designed such that the calibration sleeve can be replaced and exchanged.


Fig. 15

- ▷ To remove the calibration sleeve [3] from the tool holder [1], loosen and completely unscrew both screws [4] using the size 3 Allen key.
- ▷ Pull off the calibration sleeve [3].
- ▷ Insert a new calibration sleeve with the uncoated side onto the tool holder.

Both the calibration sleeve [3] as well as the tool holder [1] have a total of four holes.

- ▷ Screw the calibration sleeve to the tool holder only via the horizontal holes [2]. The vertical holes of the tool holder are used to affix the blank holder (see chapter 3.3.4 on page 26).

Exchanging the calibration sleeve is only necessary in case of an error in the calibration process. When exchanging the calibration sleeve, it is required to carry out the calibration process once more.

The screws of the calibration sleeve are shorter than those for the blank holder.

- ▷ Do not switch the screws or get them mixed up.

- ▷ Screw the calibration sleeve firmly to the tool holder.

## 4 Processing Ceramill Ti form blanks

### 4.1 General workflow

In order to process Ceramill Ti form blanks, a respective milling file must first be transferred to the Ceramill Motion 2 (5X).

- ▷ Switch the milling machine on.
- ▷ Start the Ceramill Motion software.
- ▷ Select the desired milling program file of the Ceramill Ti form blank and click on *Transfer File(s)*.
- ▷ Make sure that the machine is set up according to the necessary data of the program window.

If blank holder 4 is not already inserted, the *Deposit tool and change blankholder* button will flash green in the transmission window.

- ▷ Remove the blank holder from the machine (see chapter 3.2.2 on page 19).
- ▷ Clean the locating bore and the surfaces in the machine.
- ▷ Insert tooling (see chapter 4.2).
- ▷ Insert the tool holder 4 via the connection mandrel to the stop into the locating bore of the machine.


Fig. 16 Mounting the blank holder 4

- 1 Fastening screw for blank holder
- 2 Blank holder 4

- ▷ Tighten fastening screw [1] with the torque wrench until the torque wrench clicks once.
- ▷ Shut the door of the Ceramill Motion and press the left button on the machine.
- ▷ The blank loading position is advanced to.
- ▷ Mount the Ceramill Ti form blank and tighten firmly (see chapter 4.3 on page 29).
- ▷ Mount the transparent splash guard [1] to the extraction opening and lock it in place with lever [3].
- ▷ Insert the swarf sieve.
- ▷ In the Ceramill Motion software, click on *OK* to continue the transfer of the milling program.
- ▷ Start the milling process.

### 4.2 Tool placement

- ▷ Insert the Ceramill Roto TI tools with the tool shank facing outward into the tool holder of blank holder 4.
  - Place the Ceramill Roto TI 2.0 (760616) on tool location 1.
  - Place the Ceramill Roto TI 1.0 (760615) on tool location 2.
- ▷ Optionally: Insert the Ceramill Roto TI replacement tools into the tool holder of blank holder 4.
  - Place the Ceramill Roto TI 2.0 (760616) on tool location 4.
  - Place the Ceramill Roto TI 1.0 (760615) on tool location 5.

## 4.3 Blank positioning

The Ceramill blank holder for rotational milling has a mounting pin [4], onto which the Ceramill Ti form blank [8] is mounted by means of its connection geometry [6].


Fig. 17 Components for blank placement

- 1 Ceramill blank holder for rotational milling
- 2 Blank screw (already factory-mounted)
- 3 Face surface
- 4 Mounting pin
- 5 Blank nut (already factory-mounted)
- 6 Connection geometry on machine side
- 7 Millable part of blank
- 8 Ceramill Ti form blank
- 9 Labelling

### NOTE:

The blank screw [2] and the blank nut [5] can be unscrewed and removed from the hole of the mounting pin.

- ▷ When reassembling, insert the blank screw [2] into the side of the pin with the labelling [9] of the blank holder 4. This ensures that the Ceramill Ti form blank is correctly mounted (mechanical coding).


Fig. 18 Views of the Ceramill Ti form blank

- A1** Rear side of the blank!
- A2** Bottom side of the blank!
- A3** Front side of the blank!
- D1** Recess for blank nut
- D2** Recess for blank screw

The Ceramill Ti form blank has holes [D1] [D2] both on the front side [A3] as well as on the rear side [A1]. The diameter of hole [D2] is smaller than that of hole [D1]. When mounting the connection geometry [6] of the blank [8] to the mounting pin [4], the rear hole [D1] accepts the blank nut [5] and hole [D2] accepts the blank screw [2]. This ensures that the Ceramill Ti form blank cannot be mounted incorrectly, turned by 180°.

The mounting pin [4] and the face surface [3] of the Ceramill blank holder for rotational milling must be absolutely clean.

- ▷ For placement of the Ceramill Ti form blank:  
Loosen the blank screw [2] on the mounting pin [4] of the Ceramill blank holder for rotational milling.

To mount the blank [8] onto the mounting pin [4], it is **not necessary** to unscrew the blank screw [2] completely.

**NOTE:**

Incorrect assembly!

Due to the varying diameters of blank screw and blank nut, the Ceramill Ti form blank can only be mounted in one position.

► Take mounting direction into account!

► Mounting the Ceramill Ti form blank onto the mounting pin [3].


Fig. 19

**NOTE:**

Insufficient milling results!

► Tighten blank screw [2] firmly in order to prevent the Ceramill Ti form blank from becoming loose during processing.

► Firmly tighten blank screw [2] with the size 2.5 Allen key.

In the process, the blank is automatically aligned in the exact position on the Ceramill blank holder for rotational milling.

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**Amann Girrba** AG
**Amann Girrba** GmbH