Gebrauchsanweisung
Instruction Manual
Instructions d'utilisation
Istruzioni per l'uso
Instrucciones de uso

Deutsch	3 - 8
English	9 - 14
Français	15 - 20
Italiano	21 - 27
Español	28 - 34

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- Translation of the original Operating Instructions -

Table of Contents

1	Explanation of Symbols	9
2	Suitable Personnel	10
3	Product Description	10
3.1	Blank holder for glass ceramic blanks, threefold	10
3.2	Multi-block holder, twelvefold	10
4	Application	10
4.1	Loading the blank holder for glass ceramic blanks	10
4.2	Loading the multi-block holder	11
4.3	Preparing the milling machine	12
4.4	Inserting the blank holder	14
4.5	Starting the grinding program	14
5	Warranty/Exclusion of Liability	14

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a colour-backed triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- **NOTE** means that property damage can occur.
- **CAUTION** means that light to fairly serious personal injury can occur.
- **WARNING** means that serious personal injury can occur.
- **DANGER** means that serious personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 Suitable Personnel

NOTE:

The product may only be used by trained specialised personnel.

3 Product Description

The milling machine is available as four-axis model (4X) and as five-axis model (5X). These Operating Instructions apply for both models.

3.1 Blank holder for glass ceramic blanks, threefold

The blank holder for ceramic blanks is a blank holder for 4X and 5X milling machines. It is used to accommodate up to three glass ceramic blanks and the respective grinding tools.

A blank holder each is available for the 5X and 4X milling machine.

Delivery Scope

- _ Blank holder for glass ceramic blanks (threefold)
- _ Storage rack

3.2 Multi-block holder, twelvefold

The multi-block holder is a blank holder for 5X milling machines. It is used to accommodate up to 12 glass ceramic blanks and the respective grinding tools.

Delivery Scope

- _ Multi-block holder, twelvefold (5X)

4 Application

4.1 Loading the blank holder for glass ceramic blanks

The blank holder can also be inserted into the milling machine before loading it (see page 14).

- ▷ Insert the standard tools into the blank holder for glass ceramic blanks. Position the tool shank facing toward the outside.
 - Set 1.8 cutter to position 1.
 - Set 1.4 cutter to position 2.
 - Set 1.0 cutter to position 3.


Fig. 1

- ▷ Optional: Insert the sister grinding tools into the blank holder for glass ceramic blanks. Position the tool shank facing toward the outside.
 - Set 1.8 cutter to position 4.
 - Set 1.4 cutter to position 5.
 - Set 1.0 cutter to position 6.

- ▷ Depending on your selection in the CAM software, insert up to three blanks into the blank holder for glass ceramic blanks with the groove facing frontward, and secure them with the Allen key.


Fig. 2

4.2 Loading the multi-block holder

The blank holder can also be inserted into the milling machine before loading it (see page 14).

- ▷ Insert the standard tools into the multi-block holder. Position the tool shank facing toward the outside.
- Set 1.8 cutter to position 1.
 - Set 1.4 cutter to position 2.
 - Set 1.0 cutter to position 3.
 - Set 0.4 cutter to position 4.


Fig. 3

- ▷ Optional: Insert the sister grinding tools into the multi-block holder. Position the tool shank facing toward the outside.
- Set 1.8 cutter to position 6.
 - Set 1.4 cutter to position 7.
 - Set 1.0 cutter to position 8.
 - Set 0.4 cutter to position 9.
- ▷ Depending on your selection in the CAM software, insert up to 12 blanks into the multi-block holder with the groove facing frontward, and secure them with the Allen screws.


Fig. 4

4.3 Preparing the milling machine

When the milling machine is set up for dry processing:

- ▷ Convert the milling machine for wet processing (see Operating Instructions).
- ▷ Disengage lever [3] and pull off the splash guard [1] toward the front.


Fig. 5 Five-axis model


Fig. 6 Four-axis model

Legend to figure 5 and figure 6:

- 1 Splash guard
- 2 Blank holder
- 3 Lever
- A Operation button

To remove the blank holder:

- ▷ With the door closed, press operation button [A] for at least three seconds.

The current tool is placed down and the blank holder [2] returns to the change position.


Fig. 7 Five-axis model


Fig. 8 Four-axis model

Legend to figure 7 and figure 8:

- 1 Fastening screw for blank holder
- 2 Blank holder

- ▷ Loosen fastening screw [1].
- ▷ Pull off the blank holder toward the right (five-axis model) resp. toward the bottom (four-axis model).
- ▷ Clean the locating bore and surfaces.

Blank holders not in use can be placed on the storage rack and thus securely stored.


Fig. 9

4.4 Inserting the blank holder

- ▷ Check the pin of the blank holder and clean if required.
- ▷ Insert the blank holder with the pin **to the stop** in the locating bore.


Fig. 10

NOTE:

Improper handling leads to inaccurate working results or to damage of the blank-holder fixture!

- ▷ Tighten fastening screw [1] with a torque of 4 Nm.

The torque wrench included in the delivery scope of the milling machine is appropriately set.

- ▷ Tighten fastening screw [1] with with the torque wrench.

4.5 Starting the grinding program

- ▷ Briefly press operation button [A].
The machine moves to the starting position.
- ▷ Mount the splash guard [1] to the extraction opening and lock it in place with lever [3].
- ▷ Transfer the milling file for wet processing.

5 Warranty/Exclusion of Liability

Application-technical recommendations, whether given orally, in writing or in the course of practical training, are guidelines. Our products are subject to continuous further development. Subject to changes in handling and composition.

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Amann Girrbaach AG
Amann Girrbaach GmbH