# CAD/CAM Accessories

|                              |         |
|------------------------------|---------|
| <b>DE</b> Gebrauchsanweisung | 3 - 10  |
| <b>EN</b> Instruction Manual | 11 - 18 |

Kalibrierung Rotationsfräshalter Ti-Preforms (181364)  
Calibration of Blank Holder Ti-Preforms (181364)

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Table of Contents

1 Explanation of Symbols ..... 11

2 Calibration of Ti-Preforms Blank Holder ..... 12

2.1 Product description ..... 12

2.2 Delivery scope ..... 12

2.3 Required system components .... 12

2.4 Initial calibration of the Ti-Preforms blank holder ..... 13

2.4.1 Preparation of the hybrid milling machine ..... 13

2.4.2 Inserting the Ti-Preforms blank holder ..... 14

2.4.3 Measuring the measuring surfaces 15

2.4.4 Performing the final calibration milling ..... 16

2.4.5 Mounting the blank holder for rotational milling ..... 17

2.4.6 Completing the calibration process 17

2.5 Changing a calibration sleeve ..... 18

3 Spare and Wear Parts ..... 18

1 Explanation of Symbols

Warning indications

-  Warning indications in the text are marked with a triangle and boxed.
-  In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- **NOTE** means that property damage can occur.
- **CAUTION** means that light to fairly serious personal injury can occur.
- **WARNING** means that serious personal injury can occur.
- **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

| Symbol | Meaning  |
|--------|--|
| ▷      | Item of an operation description   |
| –      | Item of a list   |
| •      | Subitem of an operation description or a list                            |
| [3]    | Numbers in square brackets refer to position numbers in graphics/figures |

## 2 Calibration of Ti-Preforms Blank Holder

Before using the Ti-Preforms blank holder (181364) in the hybrid milling machine, the blank holder must be calibrated.

### NOTE:

The blank holder is individually calibrated to the respective hybrid milling machine. Afterwards, the holder can no longer be used in other hybrid milling machines.

### 2.1 Product description

Calibratable Ti-Preforms blank holder for use in hybrid milling machines. The blank holder is used for the fixation of prefabricated Ti-Preforms blanks.

The blank holder consists of two components:

- \_ Workpiece holder
- \_ Blank holder

### 2.2 Delivery scope


Fig. 1 Ti-Preforms blank holder

- 1 Workpiece holder incl. calibration sleeve
- 2 Blank holder
- 3 Additional calibration sleeve
- 4 Allen screws, M4
- 5 Allen key (size 3)
- 6 Allen key (size 2.5)

The Ti-Preforms blank holder is supplied with an already mounted calibration sleeve. A further one is enclosed.

The Ti-Preforms blank holder must be calibrated before use. For this, the calibration sleeve is face milled to match the rotation axis of the hybrid milling machine. Afterwards, the blank holder can be mounted.

### 2.3 Required system components

- \_ CAD/CAM software for the Ti-Preforms (m-plant)
- \_ Hybrid milling machine Ceramill Motion 3
- \_ Roto RFID 2.0 Ti (760672)
- \_ Roto RFID 1.0 Ti (760673)

In order to ensure the software functionality:

- ▷ Before using, install the current update of all software components.

## 2.4 Initial calibration of the Ti-Pre-forms blank holder

The calibration process consists of two steps that must be preformed one after the other.

1. Prepare the hybrid milling machine
2. Perform the calibration of the Ti-Preforms blank holder

### NOTE:

Insufficient milling results when using in other hybrid milling machines!

The blank holder is individually calibrated to the respective hybrid milling machine.

- ▷ Afterwards, do not use the holder in other hybrid milling machines.

### 2.4.1 Preparation of the hybrid milling machine

To ensure correct calibration of the Ti-Preforms blank holder:

- ▷ Carry out auto-calibration: Follow the steps in the Ceramill Motion 3 Service menu under *Calibration*.
- ▷ Perform the chuck cleaning: Follow the steps in the Ceramill Motion 3 Service menu under *Clean chuck*.

## 2.4.2 Inserting the Ti-Preforms blank holder

- ▷ Check the connection mandrel of the blank holder and clean if required.
- ▷ Insert the tool holder together with the calibration sleeve via the connection mandrel to the stop into the locating bore.


Fig. 2

- ▷ Tighten fastening screw [1] with the torque wrench until the torque wrench clicks once.


Fig. 3 Ti-Preforms blank holder

- 1 Fastening screw for tool holder
- 2 Tool holder

- ▷ Carry out a Ti-Preforms blank holder inventory run.
  - ▷ Insert tools into the tool magazine:
    - Roto RFID 2.0 Ti (760672)
    - Roto RFID 1.0 Ti (760673)
- The hybrid milling machine carries out an inventory run of the tool magazine.
- ▷ Start the hybrid milling machine software on the PC.

- ▷ Click the *Add job* button and switch to the calibration folder by clicking on the *Calibration* button.
- ▷ Select *M3 Blank Holder Calibration\_1.xnc* file and then click on the *Transfer* button to transfer data to the hybrid milling machine.


Fig. 4 Transfer "M3 Blank Holder Calibration\_1.xnc" file

- ▷ If *Start jobs after transfer* was not selected: Start job manually at the hybrid milling machine.

The calibration sleeve is now provided with 2 measuring surfaces.

The process takes approx. 20 minutes.

## 2.4.3 Measuring the measuring surfaces


Fig. 5

Two plane and opposite surfaces have been milled in the rear third of the uncoated surface of the calibration sleeve. The distance between them must be measured.

### NOTE:

Faulty calibration possible!

▷ An accurate measurement is required for a successful calibration.

- ▷ Centrally position the external jaw to one of the two surfaces and measure the distance to the opposite surface with the moving internal jaw. The milled surfaces of the calibration sleeve must be fully embraced by both caliper jaws.
- ▷ Note down the determined measuring value two digits behind the decimal point, e.g., 31.93 mm.

- ▷ If required, clean the caliper jaws.
- ▷ Set the unit of measurement on the caliper to mm.
- ▷ Slide the caliper jaws together and set to 0.

## 2.4.4 Performing the final calibration milling

▷ Start the hybrid milling machine software on the PC.

▷ Switch to calibration folder by clicking on the Calibration button.

Several program files are listed:


Fig. 6 List of program files

| MM    | FN        | MM    | FN        | MM    | FN        | MM    | FN        | MM    | FN        |
|-------|-----------|-------|-----------|-------|-----------|-------|-----------|-------|-----------|
| 32.12 | ..._32.12 | 32.07 | ..._32.07 | 32.02 | ..._32.02 | 31.97 | ..._31.98 | 31.92 | ..._31.92 |
| 32.11 | ..._32.11 | 32.06 | ..._32.06 | 32.01 | ..._32.01 | 31.96 | ..._31.96 | 31.91 | ..._31.91 |
| 32.10 | ..._32.10 | 32.05 | ..._32.05 | 32.00 | ..._32.00 | 31.95 | ..._31.95 | 31.90 | ..._31.90 |
| 32.09 | ..._32.09 | 32.04 | ..._32.04 | 31.99 | ..._31.99 | 31.94 | ..._31.94 | 31.89 | ..._31.89 |
| 32.08 | ..._32.08 | 32.03 | ..._32.03 | 31.98 | ..._31.99 | 31.93 | ..._31.93 | 31.88 | ..._31.88 |

Tab. 1 Program file - Allocation of measuring result

**MM** Measuring result in mm

**FN** File name end (M3-Holder 4-Calibration...)

▷ Select the program, whose last 4 file name digits correspond with the noted down measuring value in step 1. e.g.:

- Determined measuring value MM: 31.89 mm

- FN calibration file to be selected:

M3-Holder 4-Calibration\_31.89

▷ Transfer the respective program to the hybrid milling machine.

▷ Ensure that the machine is set up according to the necessary specifications of the program window and confirm with **OK**.

The milling data are transferred.

When *Start jobs after transfer* was not selected:

- ▷ Start job manually at the hybrid milling machine.

The surface of the calibration sleeve is now face milled.

The calibration process takes approx. 20 minutes.

## 2.4.5 Mounting the blank holder for rotational milling

### NOTE:

When assembling, the anti-rotation bolt in the blank holder must be in the provided groove in the calibration sleeve.

- ▷ Avoid jamming/wedging!

- ▷ Carefully mount the blank holder onto the calibration sleeve of the tool holder and firmly screw together with the two provided screws using a size 3 Allen key.

After screwing the blank holder and the tool holder together, a small gap will remain, which is negligible.

The blank holder for rotational milling can only be used for this machine.

- ▷ Make an appropriate mark on the blank holder.


Fig. 7

## 2.4.6 Completing the calibration process

To complete the calibration process:

- ▷ Clamp the screwed together blank holder for rotational milling into the machine.
  - Check the connection mandrel of the blank holder for rotational milling and clean if required.
  - Insert blank holder for rotational milling via the mandrel to the stop in the locating bore.
  - Tighten the fastening screw with the torque wrench.

- ▷ Mount the transparent splash guard to the extraction opening and lock it in place with the lever.

## 2.5 Changing a calibration sleeve

The blank holder for rotational milling is designed such that the calibration sleeve can be replaced and exchanged.

Exchanging the calibration sleeve is only necessary in case of an error in the calibration process. When exchanging the calibration sleeve, it is required to carry out the calibration process once more.


Fig. 8

- ▷ To remove the calibration sleeve [3] from the tool holder [1], loosen and completely unscrew both screws [4] using the size 3 Allen key.
- ▷ Pull off the calibration sleeve [3].
- ▷ Insert a new calibration sleeve with the uncoated side onto the tool holder.

Both the calibration sleeve [3] as well as the tool holder [1] have a total of four holes.

- ▷ Screw the calibration sleeve to the tool holder only via the horizontal holes [2]. The vertical holes of the tool holder are used to affix the blank holder (see chapter 2.4.5 on page 17).

The screws of the calibration sleeve are shorter than those for the blank holder.

- ▷ Do not switch the screws or get them mixed up.

- ▷ Screw the calibration sleeve firmly to the tool holder.

## 3 Spare and Wear Parts

### Art. No. Spare part

|        |                           |
|--------|---------------------------|
| 179279 | Calibration sleeve        |
| 179288 | Blank screw and blank nut |
| 760672 | Roto RFID 2,0 Ti          |
| 760673 | Roto RFID 1,0 Ti          |

Tab. 2

Manufacturer | Hersteller

**Amann Girrbach AG**  
Gewerbestrasse 10  
6841 Maeder | Austria  
Tel. +43 59 301 2100  
amanngirrbach.com

Distribution | Vertrieb D/A

**Amann Girrbach GmbH**  
Dürrenweg 40  
75177 Pforzheim | Germany  
Tel. +49 7231 957 100  
amanngirrbach.com

Made in the European Union

ISO 9001

36752\_2024-07-08