CAD/CAM Accessories

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3 - 6

EN Instruction Manual

7 - 10

Kalibrierung Werkstückhalter D-Set (181365)
Calibration of workpiece holder D-Set (181365)

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- Translation of the original Operating Instructions -

Table of Contents

1	Explanation of Symbols	7
2	Calibration of blank holder D-Set	8
2.1	Product description	8
2.2	Delivery scope	8
2.3	Required system components	8
2.4	Calibrating the blank holder	9
2.4.1	Preparation of the hybrid milling machine	9
2.4.2	Calibration of the workpiece holder	9

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

– **NOTE** means that property damage can occur.

– **CAUTION** means that light to fairly serious personal injury can occur.

– **WARNING** means that serious personal injury can occur.

– **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
–	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 Calibration of blank holder D-Set

Before using the workpiece holder D-Set (181365) in the hybrid milling machine, it must be calibrated first

NOTE:

The workpiece holder is individually calibrated to the respective hybrid milling machine. Afterwards, the holder can no longer be used in other hybrid milling machines.

2.1 Product description

Calibratable D-Set workpiece holder for use in hybrid milling machines. The blank holder is used for fixation of prefabricated D-set blanks.

The blank holder consists of two components:

- _ Blank holder
- _ Fixation plate

2.2 Delivery scope


Fig. 1 Workpiece holder D-Set Motion 3

- 1 Blank holder (uncalibrated)
- 2 Fixation plate
- 3 Allen screws M4x10 (7 pce.)

The blank holder is delivered completely assembled.

2.3 Required system components

- _ CAD/CAM software for removable prosthetics (D-Flow)
- _ Hybrid milling machine Ceramill Motion 3
- _ Roto RFID 1.2 ST (751001)

In order to ensure the software functionality:

- ▷ Before using, Install the current update of all software components.

2.4 Calibrating the blank holder

The calibration process consists of two steps that must be preformed one after the other.

1. Prepare the hybrid milling machine
2. Calibration of the blank holder

NOTE:

Insufficient milling results when using in other hybrid milling machines!

The workpiece holder is individually calibrated to the respective hybrid milling machine.

- ▷ Afterwards, do not use the holder in other hybrid milling machines.

2.4.1 Preparation of the hybrid milling machine

To ensure correct calibration of the blank holder:

- ▷ Carry out auto-calibration: Follow the steps in the Ceramill Motion 3 Service menu under *Calibration*.
- ▷ Perform the chuck cleaning: Follow the steps in the Ceramill Motion 3 Service menu under *Clean chuck*.

2.4.2 Calibration of the workpiece holder

NOTE:

Danger of tool breakage!

Calibration with the fixation plate screwed on will lead to cutter breakage and incorrect calibration.

- ▷ Loosen screws and remove fixation plate.


Fig. 2 Workpiece holder D-Set Motion 3

- ▷ Insert blank holder (without fixation plate and screws) into the material changer.

- ▷ Insert tool Roto RFID 1.2 ST (751001) into the tool magazine.

The hybrid milling machine carries out an inventory run of the tool magazine.

- ▷ Start the hybrid milling machine software on the PC.
- ▷ Switch to calibration folder by clicking on *Calibration* button.

- ▷ Select *M3_D-Set - Calibration* file and then click on the *Transfer* button to transfer data to the hybrid milling machine.


Fig. 3 Transferring the "M3_D-Set - Calibration" file

- ▷ If *Start jobs after transfer* was not selected:
Start job manually at the hybrid milling machine.
16 calibration marks are milled for fixation of the tooth blanks in the correct position.
The calibration process takes approx. 35 minutes.
- ▷ When the calibration is completed: Mount fixation plate back to blank holder with the seven screws.

The workpiece holder is now calibrated for this production machine and can immediately be used for machining D-Set blanks.

Manufacturer | Hersteller
Distribution | Vertrieb

Amann Girrbach AG
Herrschaftswiesen 1
6842 Koblach | Austria
Fon +43 5523 62333-105
Fax +43 5523 62333-5119
austria@amanngirrbach.com
www.amanngirrbach.com

Made in the European Union

ISO 9001

36753 · 2022-08-29