argovent 2

DE	Betriebsanleitung	3- 20
EN	User Manual	21- 38
FR	Mode d'emploi	39- 57
IT	Istruzioni d'uso	58- 76
ES	Modo de empleo	77- 95
PT	Instruções de operação	96- 114

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DE

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- Translation of the original Instruction Manual -

Table of Contents

1	Explanation of Symbols	22
2	General Safety Instructions	23
3	Suitable Personnel	23
4	User Manual Information	23
5	Machine Specifications	24
5.1	Delivery Scope	24
5.2	Intended Use	24
5.3	Components and interfaces	24
6	Installation	25
6.1	When switching from Argovent to Argovent 2	25
6.2	For new units with Argovent 2	28
7	Application and Operation	30
7.1	Filling the sintering bowl with pearls and powder	31
7.2	Sintering	32
7.3	Adjusting the flow of the protective gas.	33
7.4	Finalizing the sintering program	34
7.5	Results	34
7.6	Using the float sintering disc in the Argovent 2	35
8	Environmental Protection	36
9	Accessories/Spare parts	36
10	Frequently Asked Questions (FAQ)	36

1 Explanation of Symbols

Warning indications

Warning indications in the text are marked with a triangle and boxed.

In case of hazards through electricity, the exclamation mark in the warning triangle is substituted by a lightning bolt.

When there is a risk of burns from hot surfaces, the exclamation mark in the warning triangle is substituted by a heat-alert symbol.

Signal words at the beginning of a warning indication specify the type and severity of the consequences, if the measures to avert the hazard are not adhered to.

- _ **NOTE** means that property damage can occur.
- _ **CAUTION** means that light to fairly serious personal injury can occur.
- _ **WARNING** means that serious personal injury can occur.
- _ **DANGER** means that life-threatening personal injury can occur.

Important information

Important information that do not lead to hazards for humans or property damage are marked with the icon aside and are also boxed.

Other symbols in the Manual

Symbol	Meaning
▷	Item of an operation description
—	Item of a list
•	Subitem of an operation description or a list
[3]	Numbers in square brackets refer to position numbers in graphics/figures

2 General Safety Instructions

When installing, starting-up and operating the machine, always observe the following safety instructions:

WARNING:

Fire hazard!

- ▷ In case of unexpected processes in the furnace (for example, strong smoke development or odour nuisance), switch the furnace off immediately and allow the furnace to cool down to room temperature.

CAUTION:

- A defective machine can lead to malfunctions!
On detection of damage or a functional defect of the machine:
- ▷ Label the machine as defective.
 - ▷ Prevent further operation until the machine has been repaired.

NOTE:

- Increased wear of door gasket, insulation, heating rods and furnace housing!
- ▷ Avoid opening the furnace when in hot condition (above 200 °C)!

3 Suitable Personnel

NOTE:

Starting-up and operation of the machine may only be carried out by trained specialised personnel.

4 User Manual Information

NOTE:

The contents in this User Manual may possibly differ from the contents of the User Manual included with the Argotherm 2 furnace. Only the information and instructions on the use of Argotherm 2 of this User Manual are binding and to be observed. Danger of malfunctions in case of non-observance!

5 Machine Specifications

5.1 Delivery Scope

- _ Sintering compartment Argovent 2
 - Retort weight
 - Protective gas retort
 - Sintering cover
 - Separating disc
 - Sintering bowl
 - Sintering base
 - _ Floor plate of Argovent 2 (8 mm)
 - _ Gas hose with angle piece
 - _ Bow/cover tongs
 - _ Sintering pearls
 - _ Auxiliary powder
 - _ User Manual
- ▷ After unpacking, check the machine for completeness and possible transport damages. Please claim any transport damages immediately with your supplier.

5.2 Intended Use

Argovent 2 is an enclosed sintering compartment for sintering milled CoCrMo units in a protective-gas atmosphere.

This sintering compartment has been specially developed for sintering the Sintron material. Only the Sintron materials, whose properties and sintering temperatures are known, may be used. Any other or further-going use, such as the processing of other products than those intended for, as well as the handling of hazardous substances or materials hazardous to one's health, does not apply as under the intended use.

5.3 Components and interfaces


Fig. 1 Sintering compartment Argovent 2

- 1 Retort weight
- 2 Protective gas retort
- 3 Sintering cover
- 4 Separating disc
- 5 Sintering bowl
- 6 Sintering base

The Teflon hose piece connects the gas tube of the sintering base with the push-in angle piece of the gas hose.

- ▷ Do not remove the Teflon hose piece from the gas tube!

6 Installation

6.1 When switching from Argovent to Argovent 2

- ▷ Pull mains plug.
- ▷ Unscrew the six screws on the rear of the furnace and remove the back panel.


Fig. 2

- ▷ Remove gas hose at both ends.


Fig. 3

- ▷ Remove insulating stone.


Fig. 4

- ▷ Remove Argovent.
- ▷ Remove floor plate.
- ▷ Remove ceramic tube.


Fig. 5

▷ Insert floor plate:

- For Argotherm: Reinsert previously removed floor plate again.
- For Argotherm 2: Insert (8 mm) floor plate from Argovent 2 delivery scope.


Fig. 6

- ▷ Thread the hose piece of the Argovent sintering base through the hole of the rear wall and slide the sintering base to the stop into the furnace; if required, apply a light jiggling motion. The hose piece must completely project out of the rear side of the furnace.


Fig. 7

- ▷ Insert the free end of the gas hose to the internal gas connection instead of the mounted hose.


Fig. 8

- ▷ Mount the push-in angle piece of the supplied gas hose to the stop onto the hose piece of the sintering base, until the resistance of the seal has been overcome. When doing this, counterhold the sintering base in the furnace compartment.

A click noise will indicate that the push-in connection is properly engaged.


Fig. 9

If connecting the Argovent 2 fails:

- ▷ Remove the sintering base.
- ▷ Outside of the furnace, hold the sintering base by the connection pin and carefully mount the gas hose onto the hose piece of the sintering base. This slightly expands the seal and makes connecting the Argovent in the furnace easier.
- ▷ Undo the connection again, place the sintering base into the furnace and re-establish the connection.

- ▷ Mount the back panel of the furnace and screw in the six screws.


Fig. 10

- ▷ Mount the insulating stone, ensuring that the ceramic studs engage into the holes in front of the floor plate.


Fig. 11

To check if the insulating stone is properly seated:

- ▷ Carefully close the furnace lift door.
- ▷ Insert mains plug again.

6.2 For new units with Argovent 2

- ▷ Pull mains plug.
- ▷ Unscrew the six screws on the rear of the furnace and remove the back panel.


Fig. 12

- ▷ Mount rear insulating stones with the studs from the delivery scope of the furnace.
- ▷ Insert (8 mm) floor plate from Argovent 2 delivery scope.


Fig. 13

- ▷ Thread the hose piece of the Argovent sintering base through the hole of the rear wall and slide the sintering base to the stop into the furnace; if required, apply a light jiggling motion. The hose piece must completely project out of the rear side of the furnace.


Fig. 14

- ▷ Insert the free end of the supplied gas hose to the internal gas connection.


Fig. 15

- ▷ Mount the push-in angle piece of the supplied gas hose to the stop onto the hose piece of the sintering base, until the resistance of the seal has been overcome. When doing this, counterhold the sintering base in the furnace compartment.

A click noise will indicate that the push-in connection is properly engaged.


Fig. 16

If connecting the Argovent 2 fails:

- ▷ Remove the sintering base.
- ▷ Outside of the furnace, hold the sintering base by the connection pin and carefully mount the gas hose onto the hose piece of the sintering base.

This slightly expands the seal and makes connecting the Argovent in the furnace easier.

- ▷ Undo the connection again, place the sintering base into the furnace and re-establish the connection.

- ▷ Mount the back panel of the furnace and screw in the six screws.


Fig. 17

- ▷ Mount the insulating stone, ensuring that the ceramic studs engage into the holes in front of the floor plate.


Fig. 18

To check if the insulating stone is properly seated:

- ▷ Carefully close the furnace lift door.
- ▷ Insert mains plug again.

7 Application and Operation

Software update

- When switching from Argovent to Argovent 2:
Carry out software update for Argotherm 1 and 2 (import program update P1) via USB to furnace.

The software update and the corresponding User Manuals are available for download from our website.

7.1 Filling the sintering bowl with pearls and powder

NOTE:

Impairment of the sintering results!
Abrasion from the sintering bowl can selectively cause shiny areas on the sintering object.

- ▷ Do not scrape against the sintering bowl bottom with a hard object.

- ▷ Place separating disc [4] into sintering bowl [5].


Fig. 19

- ▷ Put sintering bowl on suitable precision scale and add approx. 2.5 g of auxiliary powder only on the centre of the separating disc.


Fig. 20

NOTE:

Possible reaction with sintering bowl!

- ▷ Distribute powder **only** on separating disc!
Minor quantities on the bowl bottom aside of the disc are harmless, but should possibly be avoided.

- ▷ Carefully distribute powder centrally on separating disc using a suitable instrument (e. g. metal tweezers).


Fig. 21

NOTE:

Danger of breakage of sintering bowl/cover!
With too many sintering pearls in the sintering bowl, the pressure on the sintering bowl/cover increases.

- ▷ Do not fill in more than 150 g of sintering pearls.

- ▷ Carefully pour 150 g of sintering pearls out of the can centrally onto the disc and smoothen out by gently shaking the bowl.

**NOTE:**

Possible adhering/bonding of sintering cover and sintering bowl!

When using for the first time, the powder is still very sinter-active and would bond the sintering cover and sintering bowl together.

▷ Do **not** mix in powder.

The powder may not be mixed in until after the first two sintering runs.

7.2 Sintering

▷ Place the Sintron remainder/test crown centrally onto sprinkled on pearls.

**NOTE:**

Danger of breakage of sintering bowl/cover!

When the bed of pearls is not loosened up, the pressure on the sintering bowl/cover increases.

▷ To loosen up the bed of pearls: Insert the sintering cover [3] and apply a few twisting motions until the cover abuts the floor.


Fig. 22

▷ Place sintering capsule manually or with the supplied tongs into the furnace at the centre of the sintering base.

▷ Place protective-gas retort over the sintering capsule onto the sintering base.


Fig. 23

The retort weight ensures the optimal sintering result in the Argovent 2.

▷ Because of the high weight: When using the tongs, support the retort weight manually or with additional tongs from below until placed on the retort.

▷ Position the retort weight onto the protective-gas retort.


Fig. 24

- Before starting, twist/turn the protective-gas retort on the sintering base a few times back and forth, in order to grind down possible particles/dust on the sealing surfaces.

Ensure proper seating of the retort. The sintering base must be flush with the retort.

- Close the furnace door.

NOTE:

Faulty firing run!

- Ensure that the current software is installed.

- Press button to start the furnace.

7.3 Adjusting the flow of the protective gas.

- Operate the furnace with a gas pressure of 1,5 bar.

Fig. 25 Pressure test display, Argotherm 2


Fig. 26 Argotherm manometer

Depending on the furnace type, the check values at the furnace inspection opening must be between 0.3 and 0.6 for gas value 1 [1] and above 1.0 for gas value 2 [2].


Fig. 27 Pressure test display, Argotherm 2

A leakage test at the gas outlet hole of the sintering base is not required.

7.4 Finalizing the sintering program

When the sintering program is over, a temperature of 50 °C is indicated and the compressed-air cooling is stopped.

- ▷ Switch the furnace off.
- ▷ Shut off gas cylinder at main valve.

NOTE:

Possible tilting off or falling down of the sintering bowl, due to its high weight.

- ▷ Apply bowl/cover tongs ensuring a firm and good grip.

- ▷ Using the bowl/cover tongs, carefully remove the sintering bowl out of the furnace compartment and place it down on a non-flammable surface.
- ▷ Remove the sintering cover and place it down on a non-flammable surface.
- ▷ Remove sintered frame.

When not using for longer periods, leave the Argovent components in the furnace.

7.5 Results

After the test run of the initial putting into operation, the pearls should have taken on a darker colouration, and the Sintron remainder can still have significant traces of residual oxide.

- ▷ Discard the sintering piece.

The further sintering runs should correspond with the Sintron sintering chart and not exceed the shown degree of oxidation. The results may vary within the specification range and are to be used in flawless condition. Partial discolourations within the specification range are harmless and no reason for a complaint.

The Sintron sintering chart is available for download on our website.

- ▷ In case of worse reproducible results than shown in the specification: Contact the Service Dept.

Powder exchange

For each necessary powder exchange (due to power outage, empty gas cylinder), the exchange can take place in accordance with the instructions provided with the powder. However, the powder may only be added/applied onto the separating disc (see section 7.1).

Coverings/Coatings

Whitish coverings/coatings can form between the sealing surfaces, which may facilitate the entry of air.

- ▷ To ensure tightness against leakage: After good sintering results, place retort back on base in same position.

Coverings/coatings are harmless; however, if excessive (due to heavy usage), these can be blasted off at least on the retort side using 50 µm glass pearls.

7.6 Using the float sintering disc in the Argovent 2

Sintering on pearls generally always provide a better result than float sintering with only little or no pearl/powder mixture.

NOTE:

Insufficient sintering results due to faulty disc!

- ▷ Do not use separating discs as float sintering discs, as the separating discs do not have a ground (smooth polished) surface.

For sintering large objects, the amount of pearls can be reduced in order to gain more height in the sintering bowl. It may be necessary to increase the gas pressure.

Amount of pearls	Separating disc	Gas pressure
150 g - 100 g	yes	1.5 bar
< 100 g	yes	2.0 bar
Without pearls	no	2.0 bar

Tab. 1

- ▷ Remove pearls according to desired height.
- ▷ When all pearls are removed: Remove separating disc.
- ▷ Insert float sintering disc.
- ▷ When less than 100 g of pearls are in the sintering disc: Increase gas pressure to 2.0 bar. The gas values at the inspection glass increase accordingly.


Fig. 28 Float sintering with pearls


Fig. 29 Float sintering without pearls

- 1 Float sintering disc (polished side atop)
- 2 Sintering pearls with additional powder
- 3 Separating disc
- 4 Sintering bowl

8 Environmental Protection

Packaging

In terms of packaging, the manufacturer participates in country-specific recycling systems, which ensure optimal recycling.

All packaging materials used are environmentally friendly and recyclable.

Waste equipment

Waste equipment contain materials that should be sorted for environmental-friendly recycling.

- ▷ After the service life has elapsed, dispose of the machine via the public disposal systems in an environmentally-friendly manner.

9 Accessories/Spare parts

Subject to changes.

Art. No.	Designation
178729	Auxiliary powder, Sintron
178746	Sintering base Argovent 2
178747	Sintering retort Argovent 2
178748	Sintering bowl Argovent 2
178749	Sintering cover Argovent 2
178750	Retort weight Argovent 2
178751	Separation disc Argovent 2
178752	AV2 base plate for Argotherm 2
178753	Connector Argovent 2
178754	Pliers, Argovent 2
178755	Sintering pearls Sintron Argovent 2 150 g
178756	Float sintering disc 2

Tab. 2

10 Frequently Asked Questions (FAQ)

Questions on the sintering bowl and sintering cover

Question: After the sintering run, there is breakage at the bowl and/or the cover?

Answer: It is no longer required to reduce the pearls at the bowl edge through an instrument. Twisting down the bowl into the pearls several times is definitely required, in order to loosen up the bed of pearls at the edge. Otherwise, there is acute danger of breakage of the parts! When the cover is merely pressed into the pearls (and not twisted down), the pressure can lead to damage.

Question: There is breakage/a crack in the components even though the cover was twisted properly into the pearls?

Answer: Check amount of pearls! There may never be more than the supplied 150 g of pearls filled in, in order not to additionally increase the pressure on the parts.

Question: After the first few sintering runs, the black bowls already show whitish streaks on the surface?

Answer: Whitish coverings/coatings on the SiC bowl parts are normal and can be blasted off if required (max. 50 µm glass pearls). The sintering bowl parts have been blasted off at the factory with aluminium corundum. Maintenance of the parts is not required.

Question: After use, the light sintering cover shows discolourations?

Answer: Darker discolourations (grey/green/blue) on light-coloured oxide-ceramic parts after usage are normal and harmless. After burning-in of various dusts, such as Sinteron residual powder (blue), oxide ceramics show such discolourations.

Question: After sintering, pearls adhere to the edge/rim of the sintering cover. Do these need to be completely removed each time?

Answer: Pearls adhering to the cover do not have to be completely removed. They can be left on until the next sintering run and the cover is then merely twisted into the bed of pearls again.

Question: Do the Argovent 2 components require maintenance and may they only be handled with gloves?

Answer: The Argovent 2 components can be touched with bare fingers and without gloves. Discolourations are completely harmless. The components do not require maintenance or servicing. However, it is possible to remove undesired burnt-in residue/residual oxide by blasting them off. Use a glass pearl blasting agent with a grain size of 50 µm.

Questions on discolouration/appearances on the frame

Question: What if individual selective shiny areas/spots appear on the frame?

Answer: Minor separations of the SiC bowl surface can mix-in among the pearls and adhere to the frame, which can lead to shiny spots. For this reason, covers made of white oxide ceramics were introduced, in order to protect the sintering material against complete fusion through contact reac-

tions. Selective shiny spots or areas may occur and are harmless. This appearance merely occurs on the surface and does not come from within the material. Easy removal through slight grinding or blasting off is sufficient.

Question: How can selective appearances on the frame be prevented?

Answer: Do not scrape against the bowl bottom with hard objects and afterwards allow for any abrasions to mix among the pearls.

The bowl content, however, can be cleaned:

- _ Separate the pearls from the powder,
- _ Add fresh powder
- _ Blast out bowl components before filling (Establish original condition).

Questions on sintering pearls and additional powder

Question: What about green-coloured, residual powder adhering on the bowl bottom?

Answer: Green remainders of powder on the bowl bottom or on the separating disc (after faulty sintering due to lack of argon or power loss) can be blasted off (50 µm aluminium corundum/glass pearls).

Question: What about grey flakes in the pearls?

Answer: Lightly coloured flakes in the pearls merely are sinter powder particles adhering together, which have shed from the separating disc. They can be left in the pearl mixture, mixed in or removed individually without further consideration.

Questions on the gas pressure and the gas retort

Question: May the gas pressure be increased at the gas fittings in order to achieve lower oxide results?

Answer: To facilitate oxide-free results, the gas pressure during sintering may be adjusted between 1.5 bar and 2.0 bar. A higher gas pressure leads to a higher argon consumption.

Question: What needs to be observed when assembling the components?

Answer: Pay attention to clean sealing surfaces. Due to possible residual dusts/deposits on the circumferential sealing surfaces, tolerances in the oxidation result can occur.

Questions on the float sintering disc

Question: Why does the sintering result increasingly worsen when frequent, consecutive float sintering runs are performed?

Answer: The best protective-gas coverage with the highest argon concentration is achieved when sintering with a pearls and powder mixture. For sintering runs only on the float sintering disc without pearls, the protective-gas coverage is lowest, yet sufficient. Frequent switching between float sintering and pearl sintering benefits a most uniform, oxide-free result of both procedures, and is therefore recommended.

Question: What if the separating disc would be used for float sintering instead of the float sintering disc?

Answer: The separating disc may **not** be used as float sintering disc, as it does not have a ground (smooth polished) surface.

Question: What can be done against selective burn-in spots on the surface of the float sintering disc, which can no longer be removed with normal polish?

Answer: Selective burn-in spots, for example, can be polished out using a rubber roller for zirconium oxide. This does not influence the flatness of the plate.

Question: Does the separating disc always have to be below the float sintering disc when the sintering height is not sufficient?

Answer: The separating disc with a minor amount of pearls and powder benefits the result of the float sintering. When more height is required, the float sintering disc can also be placed on the bowl bottom by itself. It separates the frame from the bowl bottom. With this procedure, a somewhat increased oxidation can occur.

General sintering tips

- _ For large-span work, never connect all crowns/pontics with connectors to the stabiliser. Connecting each second connector alternately and symmetrically is sufficient. Always connect end parts.
Recommended connector strength: rear 3.0 mm, centre 2.5 mm, front 2.0 mm. Otherwise, too many strong connectors at the frames often generate too high tensile stress, resulting in frame distortion. The crowns lift off of the float sintering disc and the bridge may warp.
- _ Normal bridges from 3 units on can also be sintered with several sintering pearls on the float sintering disc, without additionally attached stabilisers, in the case that even a deep pearl embedding is not sufficient as support (for heavy pontics next to thin-walled crowns, etc.).

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Amann Girrbach AG