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

Maintenance Instructions

1. Prevention of Surface Contamination

The milling dust of the Ceramill Sintron blanks leads to blue discolourations on the oxidation layer of the Argovent components.

Slight traces of discolouration are absolutely harmless and do not limit the functionality of the sintering structure. However, regular strong contamination through massive dust contamination leads to a burning-in effect on the oxidation layer and can occasionally lead to distortion of the Argovent components.

To prevent or minimise surface contamination through the debris and dust, the Argovent bowls should not be touched with soiled fingers, which is why the provided cotton gloves or the Argovent tongs should be used each time when putting the sintering furnace into operation.

If required, soiled components can be cleaned with the dental steam cleaner prior to the sintering process and then dried.

2. Reducing the Roughness

After several sintering cycles, the surface of the components becomes rough.

To re-establish the sliding properties of the components with respect to each other, the surfaces of the individual components need to be thoroughly sanded from time to time with suitable sanding paper.

- ① For easier handling, cut the provided sanding paper in 5cm wide stripes.

ceramill® argovent

- ② At first, thoroughly sand the interior surfaces of the protective gas retort.
- ③ The same applies for the exterior surface of the sintering bowl.
- ④ Afterwards, thoroughly sand the interior surface of the sintering bow with the sanding paper.
- ⑤ Then also sand the exterior surface of the sintering cover.
- ⑥ After treating the surfaces, thoroughly blow off the individual components or clean them with the dental steam cleaner, so that an undisturbed oxide layer can be formed again.



⚠ In no case may the surfaces be blasted off with pure aluminium corundum!

3. Maintenance Firing Run

The sanded sintering bowl must be subjected to a necessary oxidation process (without sintering pearls).

- 1) Before starting the maintenance firing run, ensure that the protective-gas cylinder is shut.
- 2) Select program 7 on the controller.
- 3) Start the program by pressing the start button.

The oxidation process takes approx. 5 hours.

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