

SOL

The perfect **SAME-DAY**
chairside solution for
DIGITAL DENTISTRY.


**COMPACT, POWERFUL,
AND ACCURATE,** ready
to print right out of the box.

Embrace SOL for your
DENTAL WORKFLOW
today, and start creating
smiles for the future.



PRODUCT FEATURES

- 
 54 LEDs calibrated for extra curing consistency
- 
 150+ Prevalidated Materials
- 
 Selective curing Extended accuracy and lifespan
- 
 49µm pixel-size Super smooth surface finish
- 
 Instant access to resin updates



PacDent
 Passion for Excellence


VOCO

Whip Mix

DMG

NextDent

DENTCA
 CAD/CAM DENTURE


ICTUS
 DENTAL


Hpdent

SHERA
 WERKSTOFF
 TECHNOLOGIE


BEGO

ARGEN

NYTE 3D

MACKLED
 RESIN EXPERTS


dentona
 www.dentona.com


SAREMCO

keystone
 industries

APPLICATIONS



Dental Model

3 per print
54 min



Temporary Crown

45 per print
51 min



Splint

3 per print
2h 07 min



Denture Base

3 per print
1h 51 min



Denture Teeth

3 per print
56 min



Surgical Guide

10 per print
53 min



IBT

6 per print
51 min



Castable Crown

45 per print
1h 25 min



Impression Tray

6 per print
3h 37 min



Gingiva

6 per print
60 min

* The printing times and pieces per print are obtained using the Large Build Platform and the optimal settings for our CURO resin range.

ACKURETTA 3D PRINTING SOLUTION

Dental clinics and labs now have the tools to effortlessly 3D print end-product restorations chairside with increased efficiency, and at a reduced cost. Ackuretta offers a validated all-in-one biocompatible solution that includes hardware, software, finishing solutions, resins, training resources, and documentation that will enable you to complete your prints in a single day.



3D Printing Resins

Ackuretta offers a complete and versatile portfolio of CURO resins specially designed for dental applications which deliver high-strength, durable prints for enhanced patient comfort.



LCD 3D Printers

Same-day chairside dental 3D printers for digital dentistry. Ackuretta's easy-to-use printers support consistent, high-quality prints for dental professionals.

Proprietary Software

ALPHA AI is the most intuitive dental-focused AI-driven 3D slicing software on the market.

ALPHA AI



Finishing Kit

Enhance your 3D printed dental applications with the Ackuretta Finishing Kit. Clean and cure your prints to ensure safe and ready-to-use outcomes that are professional-grade.

