




# Mobile Cart with Integrated Compressor



The portable oil-free compressor is suitable for any dental, veterinary, and hygiene practice.



## Features:

- Fully Integrated & Portable Mobile Cart Delivery System
- Clean Water Bottle for Self Contained System
- 3 HP Tubings for High and Low Speed Handpieces
- 120V Plug in for Compressor
- Ultra-Quiet Oilless Compressor – Low Maintenance
- 1 HP Compressor
- Great for Hygienist, Mobile Clinics, Vet Clinics
- NO Additional Plumbing Required

## COMPRESSOR

### SPECIFICATIONS

Noise Level	60db
Tank Size	2.0 Gallon Aluminum Tank (Rust Free)
Motor	1.0HP (Rated/Running -2.0HP(PEAK))
Maximum Pressure	120 psi
CFM	3.10 CFM @ 40 PSI 2.30 CFM @ 90 PSI



**QUIET & DURABLE  
OIL-FREE  
COMPRESSOR**



# Autoclave Steam Sterilization



## Model: Clave23+ - Specifications

Chamber Size	23L (9.72" x 23.6")
Number of Trays	4
Voltage	110V 50/60Hz
Overall Dimensions (WxHxD)	18.9" x 18.3" x 27.6"



## Model: Clave16+ - Specifications

Chamber Size	16L (9" x 14.2")
Number of Trays	3
Voltage	110v 50/60 Hz
Overall Dimensions (WxHxD)	17.5" x 16.1" x 23.8"

## Slow Speed Handpiece Collection

EXTERNAL WATER - NSK COMPATIBLE 



## Non Fiber Optic High Speed Handpiece



**STANDARD HEAD**  
MODEL: F-164

- ✓ >370,000 RPM
- ✓ Standard Head Ø11.2mm
- ✓ <70db
- ✓ Economical and Reliable
- ✓ Single Water Spray
- ✓ 19W Torque



21 Kenview Blvd., Unit 9, Brampton, ON, Canada L6T 5G7  
Tel. 905.799.0517 Toll Free 1.866.799.0517  
[www.flightdentalsystems.com](http://www.flightdentalsystems.com) [info@flightdentalsystems.com](mailto:info@flightdentalsystems.com)