

# CUSTOM CASSETTES

Maximize efficiency  
while protecting your  
instruments

Increase your productivity,  
saving up to six minutes  
per procedure

Heavy gauge  
stainless steel



## **Snap-In-Silicone™ Rails**

Keep your kits organized and secure  
with Snap-in-Silicone rails

- Protect sterilized instruments
- Reduce risk of sharps injuries
- Incur less breakage

U.S. Patent No. 11,324,577

# FEATURES




## Custom Cassettes


- Heavy gauge stainless steel doesn't flex or bend
- Customizable to your instrument kit
- Unique hole design for optimal flow
- Electro-polished finish for stain resistance
- Save processing time with double-sided cassettes

## Easy Crescent Glide Latch

- Fast and easy to operate, open and go
- Secure and durable



Easy to Open  
and Close



## Innovative Snap-In-Silicone™ Rails

- Uniquely designed to easily secure instruments
- Optimized instrument spacing for fast, one-handed retrieval and return
- Increase instrument exposure during sterilization with minimal contact
- Available in seven colors for easy organization

**5 YEAR UNCONDITIONAL  
GUARANTEE**

## Flexible Designs, Premium Craftsmanship

Custom Cassettes with Snap-In-Silicone™ rails can be tailored to your instrument kit. Listed below are the most popular configurations.

### Small Cassette 3" W X 8" L X 1.25" H

Available configurations:  
0-8 instruments

#### 8 INSTRUMENTS

CD1-8N1 Blue  
CD1-8N2 Teal  
CD1-8N3 Fuchsia  
CD1-8N4 Yellow  
CD1-8N5 Light Green  
CD1-8N6 White  
CD1-8N7 Purple



CD1-8N1 Shown

### Medium Cassette 5.5" W X 8" L X 1.25" H

Available configurations:  
0-16 instruments

#### 8 INSTRUMENTS

CD2-8D1 Blue  
CD2-8D1M Blue (Mat)  
CD2-8D2 Teal  
CD2-8D3 Fuchsia  
CD2-8D4 Yellow  
CD2-8D5 Light Green  
CD2-8D6 White  
CD2-8D7 Purple

#### 16 INSTRUMENTS

CD2-16N1 Blue  
CD2-16N2 Teal  
CD2-16N3 Fuchsia  
CD2-16N4 Yellow  
CD2-16N5 Light Green  
CD2-16N6 White  
CD2-16N7 Purple



### Large Cassette 8" W X 11" L X 1.25" H


Available configurations:  
0-32 instruments

#### 16 INSTRUMENTS

CD3-16D1 Blue  
CD3-16D1M Blue (Mat)  
CD3-16D2 Teal  
CD3-16D3 Fuchsia  
CD3-16D4 Yellow  
CD3-16D5 Light Green  
CD3-16D6 White  
CD3-16D7 Purple

#### 24 INSTRUMENTS

CD3-24D1 Blue  
CD3-24D1M Blue (Mat)  
CD3-24D2 Teal  
CD3-24D3 Fuchsia  
CD3-24D4 Yellow  
CD3-24D5 Light Green  
CD3-24D6 White  
CD3-24D7 Purple



Single and additional double-sided configurations are also available

Learn more at **www.nordent.com** or  
by contacting us at **800-966-7336**



# ACCESSORIES


## Snap-In-Silicone™ Rails

Available in 8 (4 top, 4 bottom) and 16 (8 top and 8 bottom). Colors available include: blue, teal, fuchsia, yellow, white, light green, and purple.




## Air Water Syringe Clips

Secure your air and water syringe tips, scissors and needle holder to the inside




## Silicone Finger Mats

Available in small, medium and large mats for larger or bulky items. Mats can be cut to your exact requirements. Only available in blue.



## Large and Small Parts Boxes

Keep small items safe and secure during sterilization.



Contact one of your **Nordent trusted advisors** today to learn how to incorporate cassettes based on your specific needs.



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 50 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM  
CRAFTSMANSHIP**



**TRUSTED  
ADVISOR**



**UNBEATABLE  
SUPPORT**

© 2023 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved. Nordent Custom Cassettes US Patent No. 11,324,577



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786

