

IMPLAMATE™ SCALERS

TITANIUM IMPLANT SCALERS

Stress-free implant
maintenance

Slender tips easily
access those
hard-to-reach areas

Safely remove
calculus and debris

Easy maintenance
without the need for
sharpening



Nordent ImplaMate™ scalers, made from the same material as implant abutments, are designed to safely glide over all implant abutment surfaces.

Solid titanium tips are thin enough to access the tightest areas, while providing the same performance as a traditional scaler.

ImplaMate™ scalers never need sharpening and can be sterilized by any method for easy maintenance.

IMPLAMATE™ TITANIUM SCALERS

PREMIUM CRAFTSMANSHIP FOR DEPENDABLE PERFORMANCE

Eliminate any worries about flexing tips and burnished calculus, with the rigidity of traditional scalers.

Implant maintenance made safe, easy and effective.

AVAILABLE PRODUCTS



CEISN67



CEIS204S



CEISBH5-6



CEISN137M



CEISLN1-2



CEISLN3-4



CEISLN5-6



A passion for premium craftsmanship and industry-leading innovations has transformed Nordent over the last 50 years into the market leader you know today. Behind the sharpest name in dentistry is a team of trusted advisors dedicated to your success and ready to provide you with unbeatable support.



**PREMIUM
CRAFTSMANSHIP**



**TRUSTED
ADVISOR**



**UNBEATABLE
SUPPORT**

© 2022 Nordent and all other marks are registered trademarks of Nordent Manufacturing Inc. All rights reserved.



T: 800.966.7336 US & Canada | T: 847.437.4780 | F: 847.437.4786
customerservice@nordent.com | www.nordent.com

