

BLZ Dental

Scanning Experience Redefine

Intraoral Scanner



INO200



INO200

Doctor's effective diagnosis and treatment assistant



Fast Scanning Speed

30 secs for full arch scan.



Higher Accuracy

Optimized algorithms for higher accuracy and precision.





Instant Heating

Power on and scan right away.



True Color Rendered 3D Model





Oral Examination Report

Dentists' smart assistant, realizing better communication with patients.



A.I. Scan



- **New Scanning Tips**

90° Scanning Tips, more scanning angles.

- **Remote Control**

Control the software remotely, achieving more hygienic and efficient scan.

- **Easy to maintain**

Detachable power cables, improving the mobility and the maintenance.



INO200

Specification

Scan Depth	23mm
Weight	270g
Size	264mm*37mm*46mm
Tip Size	20mm*17mm
Motion Control	Yes
A.I. Tech	Yes
Output Format	STL, PLY, OBJ
Power Supply	DC12V 2A
Light Source	LED
Interface	USB 3.0
Oral Examination Report	Yes
Occlusion Detect	Yes
Undercut Detect	Yes
Cervical Margin	Yes

Minimum

Computer Configuration

Operating System	Windows 10 (64-bit)
CPU	Intel Core i5 11 Gen or above
Memory	16GB
Hard Disc	512G
GPU	NVIDIA RTX 3060 6GB or above
Resolution	1920X1080 60Hz
USB Port	At least two USB 3.0 ports



BLZ Dental

BLZ Technology (Wuhan) Co., Ltd.

TEL: +86-27-6552 0255

WEB: www.blz-dental.com

Email: info@blz-dental.com

ADD: C4166, Wuhan Overseas Scholar Business Park, No.11 Dongxin Road, East Lake New Technology Development Zone, Wuhan, China



facebook



instagram