



# F20 Sintering Furnace

## Smart Operating System



### Power Failure Renewal

Continued sintering after abnormal power failure ensures efficient sintering

### Multiple Sintering Processes

Professional 20 sintering programs are available

### Large Capacity

Support 150 crowns sintering simultaneously on three sintering hold racks

## The Dental Sintering Furnace, F20



The easy entry to the denture sintering system.



Model	F20
Dimension	512X352X656mm
Weight	64Kg
Rated Voltage	220V±10% 110v±10%
Rated Frequency	50HZ
Fuse	2X~500V 32A
Rated Power	2000W
Heating Element	Four Ultrapure Silicon Molybdenum Bar
Operation Mode	7 - Inch Color Touch Screen
Maximum Temperature	1600°C
Temperature Sensor	Model-B Thermocouple

**UP3D | CONTACT US**

Web: [www.up3ds.com](http://www.up3ds.com)

Email: [info@up3d.cn](mailto:info@up3d.cn)

Tel: +86-755-26983202

