

**WEHMER**

C o r p o r a t i o n

[\*\*www.wehmer.com\*\*](http://www.wehmer.com)**PRODUCT UPGRADE NOTICE**

Dear Valued Customer,

**WEHMER** is excited to Announce the Upgrade of our Dual Model Trimmer!!

As of February 1st 2022, we'll be able to provide you with our latest DUAL MODEL TRIMMER UNIT, which comes equipped with the following:

- A 3/4hp motor. Available in 110 volts (220 volts for international use). 50/60 Hz
- Coarse Grit Diamond Wheel for faster cutting of your models -- 3/16" x 12"
- Fine Grit Diamond Wheel for a smoother finish on your models -- 3/16" x 12"
- Commercial Grade Trimming Table with knob.
- Trimming tools such as the angulation tool, angulation tool locking knob & squaring tool.
- Water spray system with factory pre-set flow control.
- Flexible 3' plastic drain hose.



<b><i>DISCONTINUED PART NUMBER</i></b>	<b><i>REPLACEMENT PART NUMBER</i></b>
108X/110	108X2/110
108X/220	108X2/220
108P/110	108X2/110
108P/220	108X2/220

The **WEHMER** Corporation appreciates and values our business relationship, and we look forward to still provide you with our Professional grade equipment & products.

Please feel free to contact us directly if you have any additional questions. Thank you for your continued support.

Best regards,

*Brenda Mendoza*

Brenda Mendoza

Finance & Administrative Director

**1151 N. Main Street Lombard, IL 60148 USA**



**Toll Free: 1-800-323-0229 | Direct: 630-424-1877**