




























VALIDATED DEVICES | COMPATIBILITY

Pac-Dent is dedicated to continuing expansion of the list of devices that are compatible with our Rodin™ 3D Resins. We believe in giving you freedom of choice when it comes to your preferred printing device. See below for our most up-to-date validation and verification status chart.

Validated Status		✓ Verified Compatibility			⌚ Scheduled Validation		● Incompatible	
Compatible Printers		 Rodin® Sculpture 2.0	 Rodin® Titan	 Rodin® Sculpture 1.0	 Rodin® Denture Base 1.0	 Rodin® Denture Base 2.0	 Rodin® Splint 2.0	 Rodin® Clear Bite Splint
	Ackuretta SOL	✓	✓	✓	✓	✓	✓	✓
	Asiga® MAX (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Asiga® Pico 2 (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Asiga® Pro 4K (385 nm)	✓	✓	✓	✓	✓	✓	✓
	Carbon® 3D (M1, M2, M3)	⌚	⌚	✓	⌚	⌚	⌚	⌚
	DentaFAB SEGA	✓	✓	✓	✓	⌚	⌚	✓
	Flashforge® Focus 6K	●	✓	✓	✓	✓	✓	✓
	LuxCreo iLux Pro Dental	✓	✓	⌚	⌚	✓	⌚	⌚
	Nexa Xip	⌚	⌚	✓	⌚	✓	✓	⌚
	Phrozen® Sonic Mini 4K	✓	✓	✓	●	✓	✓	●
	Phrozen® Sonic Mini 8K	✓	✓	✓	⌚	✓	✓	✓
	Phrozen® Sonic Mini 8K S	✓	✓	✓	⌚	✓	✓	⌚
	Phrozen® Sonic 4K	✓	✓	✓	●	✓	✓	●
	Phrozen® Sonic 4K XL	⌚	⌚	✓	●	✓	✓	●
	RayShape Shape 1+Dental	✓	✓	✓	✓	✓	✓	✓
	RayShape Edge E1	✓	✓	✓	✓	✓	✓	✓
	Shining 3D L4D/L4K	✓	✓	✓	✓	✓	✓	✓
	Shining AccuFab-CEL	⌚	✓	✓	✓	✓	✓	✓
	Whip Mix VeriEKO	⌚	⌚	✓	✓	⌚	⌚	⌚















VALIDATED DEVICES | COMPATIBILITY







Pac-Dent is dedicated to continuing expansion of the list of devices that are compatible with our Rodin™ 3D Resins. We believe in giving you freedom of choice when it comes to your preferred printing device. See below for our most up-to-date validation and verification status chart.

Validated Status		✔ Verified Compatibility			⌚ Scheduled Validation		● Incompatible		
Compatible Printers		 Rodin® Flexisplint	 Rodin® Custom Bite Tray	 Rodin® Surgical Guide 2.0	 Rodin® Surgical Guide 1.0	 Rodin® Try-in Denture	 Rodin® Ortho IBT	 Rodin® Model	 Rodin® Castable
	Ackuretta SOL	⌚	✔	✔	✔	✔	✔	✔	●
	Asiga® MAX (385 nm)	✔	✔	✔	✔	✔	✔	✔	✔
	Asiga® Pico 2 (385 nm)	✔	✔	✔	✔	✔	✔	✔	✔
	Asiga® Pro 4K (385 nm)	✔	✔	✔	✔	✔	✔	✔	✔
	Carbon® 3D (M1, M2, M3)	⌚	✔	⌚	⌚	⌚	⌚	⌚	⌚
	DentaFAB SEGA	⌚	✔	✔	⌚	✔	✔	✔	✔
	Flashforge® Focus 6K	⌚	✔	✔	✔	✔	✔	✔	●
	LuxCreo iLux Pro Dental	⌚	⌚	⌚	⌚	✔	⌚	⌚	⌚
	Nexa Xip	⌚	⌚	✔	⌚	⌚	⌚	✔	⌚
	Phrozen® Sonic Mini 4K	⌚	✔	✔	●	✔	✔	✔	●
	Phrozen® Sonic Mini 8K	⌚	⌚	✔	⌚	✔	✔	✔	⌚
	Phrozen® Sonic Mini 8K S	⌚	⌚	✔	⌚	✔	⌚	✔	⌚
	Phrozen® Sonic 4K	⌚	✔	✔	●	✔	✔	✔	●
	Phrozen® Sonic 4K XL	⌚	✔	✔	●	✔	✔	✔	●
	RayShape Shape 1+Dental	⌚	✔	✔	✔	✔	✔	✔	⌚
	RayShape Edge E1	⌚	✔	✔	✔	✔	✔	✔	⌚
	Shining 3D L4D/L4K	⌚	✔	✔	✔	✔	✔	✔	✔
	Shining 3D AccuFab-CEL	⌚	✔	✔	✔	✔	✔	✔	✔
	Whip Mix VeriEKO	⌚	⌚	⌚	⌚	⌚	⌚	⌚	⌚

COMPATIBLE LIGHT CURING DEVICES

Pac-Dent is dedicated to continuing expansion of the list of devices that are compatible with our Rodin™ 3D Resins. We believe in giving you freedom of choice when it comes to your preferred printing device. See below for our most up-to-date validation and verification status chart.

Validated Status		✓ Verified Compatibility			⌚ Scheduled Validation		● Incompatible	
Compatible Light Curing Devices								
	Rodin® Sculpture 2.0	Rodin® Titan	Rodin® Sculpture 1.0	Rodin® Denture Base 1.0	Rodin® Denture Base 2.0	Rodin® Splint 2.0	Rodin® Clear Bite Splint	Rodin® Flexisplint
 Phrozen Curing Station	⌚	⌚	⌚	⌚	⌚	⌚	⌚	⌚
 Dreve PCU LED N2	✓	✓	✓	✓	✓	✓	⌚	⌚
 Ackuretta Curie Plus	✓	✓	✓	✓	✓	✓	⌚	⌚
 Ackuretta Curie	✓	✓	✓	✓	✓	✓	⌚	⌚
 Formlabs® Form Cure	✓	✓	✓	✓	✓	✓	⌚	⌚
 Otoflash	✓	✓	✓	✓	✓	✓	✓	✓

Compatible Light Curing Devices	Rodin® Custom Bite Tray	Rodin® Surgical Guide 2.0	Rodin® Surgical Guide 1.0	Rodin® Try-in Denture	Rodin® Ortho IBT	Rodin® Model	Rodin® Castable	Rodin® Palette
 Phrozen Curing Station	⌚	✓	⌚	⌚	⌚	⌚	⌚	⌚
 Dreve PCU LED N2	✓	✓	✓	✓	✓	✓	✓	⌚
 Ackuretta Curie Plus	✓	✓	✓	✓	✓	✓	✓	⌚
 Ackuretta Curie	✓	✓	✓	✓	✓	✓	✓	⌚
 Formlabs® Form Cure	✓	✓	✓	✓	✓	✓	✓	⌚
 Otoflash	✓	✓	✓	✓	✓	✓	✓	✓
FOR BEST AESTHETICS AND MECHANICAL PROPERTIES. MOST RECOMMENDED								