

On-X Teeth UV 1.0

- ✓ High flexural strength
(Outstanding occlusal load capacity and fracture resistance)
- ✓ High Abrasion Resistance
- ✓ Anti-corrosion



All-on-X
(Hybrid denture)



Specifications 500 g/Bottle

Machine Compatibility



A2D



ChairSide Pro

Applications

Properties	Standard	Units	Heygears on-x teeth	Competitors
Flexural Strength	ASTM D790	MPa	170	/
Flexural Strength	ISO20795	MPa	117	72
Flexural Modulus	ISO20795	MPa	3200	2388
Work of Fracture	ISO20795	J/m ²	2636	1033
Determination of MMA Content	ISO20795	/	Pass	Pass
Hardness	ASTM D2240	Shore D	91	90