



**Bellagio
Operatory Package
#BP5015-L700LED**
(Retail: \$15,170.00)



Standard Features

Mirage Hydraulic Patient Chair

- Cast aluminum base
- Chair top swivel w lock
- Foot control
- 2 Presets
- Double articulating headrest


Bellagio Delivery System

- Built in piezo scaler / LED curing light
- Integrated touch pad
- 3 hp automatic block
- Integrated timer
- Hve, Se, 3 way syringe


Vivid LED Light

- 35,000 LUX
- Manual toggle ON/OFF
- Motion sensor ON/OFF
- Dimmer function
- Composite cure mode

Standard Vinyl Upholstery Color Selection




Morbern Allante FFree Standard Vinyl



AF colors are subject to upcharge for Laguna Upholstery

Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade
Retail: \$490

Doctor Stool Upgrade
Retail: \$148

Assistant Stool Upgrade
Retail: \$170



**Mirage 2.0
Operatory Package
#MP2015-L600LED-2.0**
(Retail: \$13,080.00)



Standard Features

Mirage Hydraulic Patient Chair

- Steel cut base
- Chair top swivel w lock
- Foot control
- 2 Presets
- Double articulating headrest


Mirage Delivery System

- 3 hp automatic block
- Touch Pad controls
- 3 hp automatic block
- Hve, Se, 3 way syringe
- Water bottle system

Lucent LED Light

- 35,000 LUX
- Manual toggle ON/OFF
- Removable handles
- Dimmer function


Standard Vinyl Upholstery Color Selection



Morbern Allante FRee Standard Vinyl



Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade
Retail: \$490

Doctor Stool Upgrade
Retail: \$148

Assistant Stool Upgrade
Retail: \$170





Mirage 2.0
Swing Mount Package
#MSP3500-LED550-2.0
(Retail: \$14,256.00)



Standard Features

Mirage Hydraulic Patient Chair

- Steel cut base
- Chair top swivel w lock
- Foot control
- 2 Presets
- Double articulating headrest


Mirage Swing Mount Delivery System

- 3 hp automatic block
- Touch pad control
- 3 hp automatic block
- Hve, Se, 3 way syringe
- Water bottle system

Lustrous LED Light

- 30,000 LUX
- Manual toggle ON/OFF
- Motion Sensor ON/OFF
- Dimmer function
- Removable handles


Standard Vinyl Upholstery Color Selection



Morbern Allante FFree Standard Vinyl



Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade
Retail: \$490

Doctor Stool Upgrade
Retail: \$148

Assistant Stool Upgrade
Retail: \$170





**Laguna 2.0
Operatory Package
#LP2015-L550LED-2.0**
(Retail: \$10,080.00)



Standard Features

Laguna Electromechanical Patient Chair

- Steel Cut base
- Chair safety switch
- Foot control
- 2 Presets
- Double articulating headrest

Mirage Delivery System

- 3 hp automatic block
- Touch pad controls
- 3 hp automatic block
- Hve, Se, 3 way syringe
- Water bottle system

Lustrous LED Light

- 30,000 LUX
- Manual toggle ON/OFF
- Motion sensor ON/OFF
- Removable handles
- Composite cure mode

Standard Vinyl Upholstery Color Selection



Morbern Allante FRee Standard Vinyl



Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade
Retail: \$490

Doctor Stool Upgrade
Retail: \$148

Assistant Stool Upgrade
Retail: \$170



**Laguna 2.0
Orthodontic Package
#LOP2000-L550LED-2.0**
(Retail: \$9,416.00)



Standard Features

Laguna Electromechanical Patient Chair

- Steel cut base
- Chair safety switch
- Foot control
- 2 Presets
- Double articulating headrest


Mirage Orthodontic System

- 3 hp automatic block
- 5 Locking drawers
- Utility center w/regulators
- Hve, Se, 3 way syringe
- Pull out writing surface


Lustrous LED Light

- 30,000 LUX
- Manual toggle ON/OFF
- Motion sensor ON/OFF
- Removable handles
- Composite cure mode


Standard Vinyl Upholstery Color Selection



Morbern Allante FRee Standard Vinyl



Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade
Retail: \$490

Doctor Stool Upgrade
Retail: \$148

Assistant Stool Upgrade
Retail: \$170



#4000-2.0-CB
Mirage 2.0 Hydraulic Chair
Cast Base
Retail: \$8,590

#3000-2.0-CB
Mirage 2.0 Hydraulic Chair
Cast Base (No Armrest)
Retail: \$8,390



#4000-2.0
Mirage 2.0 Hydraulic Chair
Retail: \$7,390

#3000-2.0
Mirage 2.0 Hydraulic Chair
(No Armrest)
Retail: \$7,190



#L2000-2.0
Laguna 2.0 Electromechanical
Patient Chair
Retail: \$4,590

#L3000-2.0
Laguna 2.0 Electromechanical
Patient Chair (no armrests)
Retail: \$4,390



#LOP2000-L550LED-2.0
Laguna 2.0 Orthodontic Package

Retail: \$9,416



#PC2630
Portable Dental Unit

- 2 hp manual selection
- 3 way syringe
- Water bottle system
- Waste container
- 3/4 hp compressor
- Air venturi vacuum

Retail: \$2,390



#PC2930
Portable Dental Unit

- 2 hp automatic selection
- **LED 39 curing light**
- **Built in piezo scaler**
- 1 year warranty

Retail: \$3,300



#PC2750
Tripod LED Light

- LED illumination
- Tripod with carry bag
- Intensity adjustment
- 2 year warranty

Retail: \$378



#PC2700
Portable Dental Chair Package
Retail: \$1,996

Model	Retail
• PC 2720 Patient Chair	\$1,090
• PC 2740 Dr. Stool	\$480
• PC 2760 LED Light	\$260
• PC 2780 Tray / Spittoon	\$290
• PC 2790 Carry Bag	\$190
• PC 2795 Stool Carry Bad	\$100

#L7601-LED

Luminous LED Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle switch on / off
- 2 year warranty

Retail: \$2,900




#L7605-LED

Luminous LED Light

- **Sensor on / off**

Retail: \$3,190



#L7901-LED

Luminous LED Ceiling mount

Retail: \$4,390

#R6101-LED

Radiant LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle switch on / off
- 2 year warranty

Retail: \$1,990



#R6105-LED

Radiant LED Operatory Light

- **Sensor on / off**

Retail: \$2,190

#L600-LED

Lucent LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- Toggle on / off
- 2 year warranty

Retail: \$1,500



#L690-LED

Lucent LED Operatory Light

- **Ceiling mounted**

Retail: \$2,160



#L700-LED

Vivid LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- **Sensor +Toggle on / off**
- Composite cure mode
- 2 year warranty

Retail: \$1,590

#L795-LED

Vivid LED Operatory Light

- **Ceiling mounted**

Retail: \$2,250

#L625-LED

Lucent LED Operatory Light

- Post mounted
- LED illumination
- Lux: 3,000 - 35,000
- Intensity adjustment
- **Sensor on / off**
- 2 year warranty

Retail: \$1,700

#L695-LED

Lucent LED Operatory Light

- **Ceiling mounted**

Retail: \$2,360

#L550-LED

Lustrous LED Operatory Light

- Post mounted
- LED illumination
- Lux: 8,000 - 30,000
- Intensity adjustment
- Toggle / sensor
- 2 year warranty

Retail: \$1,300

#L570-LED

Lustrous LED Operatory Light

- **Ceiling mounted**

Retail: \$1,960



#5000-2.0

Bellagio Delivery System

- 3 hp position w Touch pad
- Cuspidor
- Assistants Package w Touch pad
- Hve / SE
- Wet dry foot control
- **Built in Piezo scaler**
- **Built in LED curing light**
- 3 way syringe
- Water bottle system
- 5 year warranty

Retail: \$5,690



#2015-2.0

Mirage 2.0 Delivery System

- 3 hp position w touch pad
- Assistants Package w touch pad
- HVE / SE
- Wet dry foot control
- 3 way syringe
- Water bottle system
- 5 year warranty

Retail: \$4,190



#5015-2.0

Bellagio Delivery System (no cuspidor)

Retail: \$4,990

#2000-2.0

Mirage 2.0 Delivery System (with cuspidor)

Retail: \$4,890

#MC500-2.0

Mirage Mobile Cart

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center with j-box
- 5 year warranty

Retail: \$2,500



#MR4500-2.0

Mirage 2.0 Rear Delivery with Vacuum Package

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty
- **HVE**
- **SE**
- **Vacuum trap**

Retail: \$2,990



#MC501-2.0

Mirage Mobile Cart

- 3 hp position
- 3 way syringe, HVE + SE
- Wet dry foot control
- Water bottle system
- 5 year warranty

Retail: \$2,800

#MS3000

Mirage Side Delivery Cabinet Mount

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty



Retail: \$3,300

#MR4000-2.0

Mirage 2.0 Rear Delivery

- 3 hp position
- 3 way syringe
- Wet dry foot control
- Water bottle system
- Utility center
- 5 year warranty

Retail: \$2,200




#OMC2375 Orthodontic Delivery

- Three way syringe
- HVE / SE valves
- Wet / dry foot control
- Water bottle system
- Utility center junction box
- 5 Storage drawers
- 1 year warranty
- **Requires central air and vacuum**



Retail: \$3,190



OMC-ST81



OMC-ST82



OMC-ST85



OMC-ST86

Retail: \$420

#OMC2375CV Orthodontic Delivery

- Three way syringe
- HVE / SE valves
- Wet / dry foot control
- Water bottle system
- Waste bottle system w/ solids collector
- 5 Storage drawers
- 3/4 hp oil-less motor
- 1 year warranty
- **Self contained**



Retail: \$4,700

#OMC2375CV-SL Orthodontic Delivery

- **Built in Piezo scaler**
- **Built in LED 39 curing light**

Retail: \$5,600



#MC-3600CV

Mirage Mobile Delivery Cart

- Self contained
- Three way syringe
- Wet / dry foot control
- Water bottle system
- Waste bottle system
- Hve / SE valves
- 1 year warranty



Retail: \$4,270

#MC-3600CV-SL

Mirage Mobile Delivery Cart

- Self contained
- Three way syringe
- Wet / dry foot control
- Water bottle system
- Waste bottle system
- Hve / SE valves
- **Built in Piezo scaler**
- **Built in LED 39 curing light**
- 1 year warranty



Retail: \$5,170

#PC2635

Portable Delivery System

- Three way syringe
- Water bottle system
- 2 HP positions
- Disc style foot control
- Carry handle
- Pull handle
- 1 year warranty



Retail: \$1,790

#PC2530

Mobile Vacuum System

- 2 HVE valves
- 1 SE valves
- Solids collector trap 2 1/8"
- 7/8 hp vacuum motor
- 2L Waste container
- 4 large durable casters
- 1 year warranty



Retail: \$2,100



- #DR-7000
Solace Doctor's Stool**
- 140 mm gas cylinder
 - Backrest tilt adjustment
 - Wide contoured seat pan
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$818



- #DR-1102
Mirage Doctor's Stool**
- 140mm gas cylinder
 - Backrest tilt adjustment
 - Wide seat pan
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$650



- #HY-5107
Doctor's Stool**
- 140 mm gas cylinder
 - Backrest tilt adjustment
 - Triangle seat pan
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$690



- #SS-5105
Mirage Saddle Stool**
- 200 mm gas cylinder
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$690



- #DR-5045
Laguna Doctor's Stool**
- 140mm gas cylinder
 - Backrest tilt adjustment
 - Wide seat pan
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$480



- #AS-1101
Mirage Assistant's Stool**
- 200 mm gas cylinder
 - Backrest tilt adjustment
 - Wide seat pan
 - Body support bar
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$818



- #AS-5047
Laguna Assistant's Stool**
- 200 mm gas cylinder
 - Body support bar
 - Circular seat pan
 - 5 dual wheel casters
 - 1 year warranty

Retail: \$590




- #XR-6101
X-Ray Chair**
- Adjustable headrest
 - Pneumatic pump for base
 - Arm rest
 - Leg rest
 - 1 year warranty

Retail: \$1,550


Standard Vinyl Upholstery Color Selection



Morbern Allante FFree Standard Vinyl



Optional NaguaSoft Vinyl Upholstery Upgrade:



Dental Chair Upgrade

Retail: \$490

Doctor Stool Upgrade

Retail: \$148

Assistant Stool Upgrade

Retail: \$170

Mechanical Room Combo Bundle



DC8228



WV3

SuperbAIR #DC8124 2HP Single Air Compressor (3 User)

OR

SuperbAIR #DC8228 4HP (Twin) Air Compressor (7 User)

+

Any model Superb VAC Vacuum Pump



DC8124



WV10

Mechanical Room

TPC Superb Vac Wet Ring Vacuum System

Model	Retail
WV3 (Single) 1 HP 110V or 220V	\$3,190
WV3R (Single) 1 HP 220V w/ Recycler	\$3,590
WV5 (Single) 2 HP 220V	\$3,590
WV5R (Single) 2 HP 220V w/ Recycler	\$3,990
WV6 (Twin) 1 HP 220V	\$6,590
WV6R (Twin) 1 HP 220V w/ Recycler	\$6,990
WV10 (Single) 2 HP 220V	\$7,390
WV10R (Single) 2 HP 220V w/ Recycler	\$7,790



TPC Superb Air Oil-Less Air Compressor Model

Model	Retail
DC8112 (Single), 1 HP 110V (1 Users)	\$5,950
DC8124 (Single), 2 HP 220V (3 Users)	\$6,590
DC8224 (Twin) , 2 HP 220V (3 Users)	\$7,990
DC8226 (Twin) , 3 HP 220V (5 Users)	\$9,590
DC8228 (Twin) , 4 HP 220V (7 Users)	\$11,590



Mobile Cabinet Systems



#TMC-120 Assistant's North Carolina Mobile Cabinet
 • 21.5"W x 19"D x 32"H
 • Color: White or light grey
 • Solid top with pull handle
 • Amalgamator well with drop leaf door
 • 3 instrument drawers
 • 1 storage drawer
 • 1 year warranty
 Shipping dimensions: 25" x 23" x 36" 112 lbs
Retail Price: \$990



#TMC-140 Doctor's Mobile Cabinet
 • 21.5"W x 19"D x 32"H
 • Color: White or light Grey
 • Solid top with pull handle
 • 4 instrument drawers
 • 3 storage drawer
 • 1 year warranty
 Shipping dimensions: 27" x 25" x 36" 116 lbs
Retail Price: \$990



#TMC-160 Assistant's Mobile Cabinet
 • 21.5"W x 19"D x 32"H
 • Color: White or light Grey
 • Solid top with pull handle
 • Sliding counter top with amalgamator well
 • 1 instrument drawers
 • 2 storage drawer
 • 1 year warranty
 Shipping dimensions: 25" x 23" x 35" 115 lbs
Retail Price: \$990



#TMC-180 Utility Mobile Cabinet
 • 14.5"W x 16"D x 32"H
 • Color: White or light Grey
 • 1 instrument drawers
 • 2 Adjustable shelves
 • 1 year warranty
 Shipping dimensions: 24" x 24" x 37" 72 lbs
Retail Price: \$478

