

VINTAGE ZR ■ VINTAGE ZR Uni-Layer



Porcelain system for all available zirconia frameworks

Vintage ZR is a fine structured, leucite reinforced microporcelain for metal-free restorations. The natural feldspar-based structure guarantees a secure bond to all zirconia frameworks and a very high stability and strength of the aesthetic restorations.



Due to the light optical properties of Vintage ZR, the porcelain restorations appear natural and expressive.

BENEFITS

- High-fusing Opaque Liner (960 °C) secure the bond to the zirconia framework and are the shade basis of the restoration
- Impressive aesthetics by natural translucency and opalescence
- Excellent shade match to natural teeth and to the VITA* Classical Shade System
- Extensive range of translucent incisal and effect powders
- Low investment cost due to modular system components
- Constant CTE of about $9.3 \times 10^{-6} \text{K}^{-1}$, even after several firings

VINTAGE ZR Uni-Layer


Fully anatomical anterior and posterior restorations can be built up with only one layer with this innovative porcelain.

Perfectly balancing light diffusion and chroma, this porcelain combines the optical characteristics of dentine and incisal powders in one material – for a harmonious shade progression and the illusion of depth.

- Quick and easy single-layer technique
- Simple correction options
- Easy shade matching with Vintage Art porcelain stains
- Can be combined with the opalescent and translucent materials of the Vintage ZR system depending on the demand and situation

Indications

Facing of all crown and bridge constructions made of zirconia



The restoration with Vintage ZR Uni-Layer can be produced in higher or lower translucency

VINTAGE ZR - Sets



AB Set

AB Set

Components for aesthetic all-ceramic restorations in the shades A1-A4, RootA and B1-B4

- 10x Opaque Liner, 5 g each • 10x Body, 15 g each
- 2x Cervical, 15 g each • 5x Opaque Dentin, 15 g each
- 4x Opal Enamel, 15 g each • Opal Effect, 15 g
- 2x correction porcelain (Add-On B/Add-On T), 15 g each
- Vintage AL/ZR Opaque Liner Liquid, 3 ml • Vintage Modelling Liquid, 50 ml
- colour table • instructions for use



CD Set

CD Set

Practical addition in the shades C1-C4 and D2-D4 to complete the AB Set

- 7x Opaque Liner, 5 g each • 7x Body, 15 g each • 2x Cervical, 15 g each
- 3x Opaque Dentin, 15 g each • Vintage AL/ZR Opaque Liner Liquid, 3 ml
- colour table • instructions for use

Enamel Effect Set

Opalescent and translucent incisal powders for individual effects in the incisal area

- 8x Opal Effect, 15 g each • 5x Enamel Effect, 15 g each
- 5x Cervical Translucent, 15 g each • colour table
- instructions for use

Color Effect Set

For individual effects in dentine and cervical areas

- 3x Opaque Liner Effect, 5 g each • 8x Color Effect, 15 g each
- 2x Gum, 15 g each • colour table • instructions for use

Whitening Set

Additional shades for brighter restorations

- 3x Opaque Liner, 5 g each • 3x Body, 15 g each
- 2x Opaque Dentin, 15 g each • 2x Opal Enamel, 15 g each
- Opal Effect, 15 g • Vintage AL/ZR Opaque Liner Liquid, 3 ml
- colour table • instructions for use

PN	Quantity	Product
9000	1 pc.	Vintage ZR AB Set
9001	1 pc.	Vintage ZR CD Set
9002	1 pc.	Vintage ZR Whitening Set
9003	1 pc.	Vintage ZR Enamel Effect Set
9005	1 pc.	Vintage ZR Color Effect Set

VINTAGE ZR - Single Components

Shade	PN 5 g
OPAQUE LINER	
W10	9030
W20	9031
W30	9032
A10	9011
A20	9012
A30	9013
A3.50	9014
A40	9015
RootAO	9016
B10	9017
B20	9018
B30	9019
B40	9020
C10	9021
C20	9022
C30	9023
C40	9024
D20	9026
D30	9027
D40	9028
OPAQUE LINER EFFECT	
OM-Y	9034
OM-LP	9035
OM-DP	9036

Shade	PN 15 g	PN 50 g
MARGIN		
W1M	9218	
NM	9237	
A3M	9220	
RootAM	9222	
B2M	9224	
B4M	9226	
C2M	9229	
C4M	9231	
D3M	9233	
MARGIN EFFECT		
CLM	9236	
LPM	9238	
CERVICAL		
AC	9173	9179
BC	9174	9180
CC	9175	9181
DC	9176	9182
BODY		
W1B	9097	9122
W2B	9098	9123
W3B	9099	9124
A1B	9077	9102
A2B	9078	9103
A3B	9079	9104
A3.5B	9080	9105
A4B	9081	9106
RootAB	9082	9107
B1B	9083	9108
B2B	9084	9109
B3B	9085	9110
B4B	9086	9111
C1B	9088	9113
C2B	9089	9114
C3B	9090	9115
C4B	9091	9116
D2B	9093	9118
D3B	9094	9119
D4B	9095	9120

Shade	PN 15 g	PN 50 g
OPAQUE DENTIN		
OD-W1	9041	9060
OD-N	9039	9058
OD-A3	9043	9062
OD-RootA	9045	9064
OD-B2	9047	9066
OD-B4	9049	9068
OD-C2	9051	9070
OD-C4	9053	9072
OD-D3	9055	9074
COLOR EFFECT		
MP	9199	
MY	9200	
Mlv	9201	
R	9202	
Y	9203	
O	9204	
G	9205	
W	9206	
OPAL ENAMEL		
Opal 56	9127	9143
Opal 57	9128	9144
Opal 58	9129	9145
Opal 59	9130	9146
Opal 60	9131	9147
OPAL EFFECT		
Opal T	9133	9149
Opal SL	9134	9150
Opal WE	9135	9151
Opal MI	9136	9152
Opal OC	9137	9153
Opal AM-R	9138	9154
Opal AM-Y	9139	9155
Opal AM-V	9140	9156

Shade	PN 15 g	PN 50 g
ENAMEL EFFECT		
BT	9159	9166
OT	9160	9167
PT	9161	9168
GT	9162	9169
T-Glass	9163	9170
CERVICAL TRANSLUCENT		
CT-CL	9185	9192
CT-W	9186	9193
CT-A	9187	9194
CT-B	9188	9195
CT-R	9189	9196
GUM		
Gum-R	9207	
Gum-P	9208	
Gum-LP	9209	
Gum-DP	9210	
Gum-V	9211	
Gum-Or	9212	
CORRECTION		
Add-On B	9213	
Add-On T	9214	
CPM Fine	9215	



Uni-Layer Set

Contains all components for aesthetic restorations in the shades A2 and A3

- Opaque Liner A20, 5 g • Opaque Liner A30, 5 g
- Vintage ZR Uni-Layer HT, 50 g • Vintage ZR Uni-Layer LT, 50 g
- Vintage AL/ZR Opaque Liner Liquid, 3 ml
- Vintage Mixing Liquid HC, 50 ml • step card • instructions for use

PN	Quantity	Product
9247	1 pc.	Vintage ZR Uni-Layer Set (A2 / A3)
9241	50 g	Vintage ZR Uni-Layer HT
9244	50 g	Vintage ZR Uni-Layer LT

VINTAGE ZR - Liquids

PN	Quantity	Product
8011	3 ml	Vintage AL/ZR Opaque Liner Liquid
9611	50 ml	Vintage Mixing Liquid
9612	500 ml	Vintage Mixing Liquid
9613	50 ml	Vintage Mixing Liquid HC
9614	500 ml	Vintage Mixing Liquid HC
9615	50 ml	Vintage Margin Hardening Liquid
7642	3 ml	Vintage CPM Modelling Liquid
8012	7 ml	Vintage Margin Porcelain Isolation Liquid (pen)