

Ray Foster Dry Model Trimmer



*Catalog MTD1
Dry Model Trimmer*

Features:

Rigid cutting surface
gives a flat, smooth finish

Simplified mechanics
(no noisy pulleys or belts) results in a
smooth running, quiet machine

Uses fast cutting PSA abrasive discs

Heavy-duty, capacitor-start,
1/3 HP motor

Adjustable work table

Six foot hose and four additional
abrasive discs included

Self contained model trimming system
when connected to the Ray Foster
Cyclone Dust Collector

Proven Ray Foster quality
and durability

*Catalog MTD1C
Dry Model Trimmer
with Cyclone Dust
Collector*



Specifications

14" high, 14" wide, 15" deep
115 volt, 60 hz standard
230 volt, 50/60 hz optional for export
Wheel size 10" diameter
1-1/2" O.D. suction connection
Net weight 32 lbs
Shipping weight 40 lbs

Produce high quality,
precise models with-
out water.

When dry trimming
of stone models is
required, as with refrac-
tory investments, or
preferred, the Ray Foster
Dry Model Trimmer is
the right choice.

www.fosterdental.com