



UP3D Desktop Scanner System Requirements

Minimum System Requirements:

CPU: Intel i5-2500 or better (4 CPU cores or more.)

Graphics: Nvidia 2 GB Dedicated Video RAM

RAM: 8GB DDR4

Storage: 240GB SSD (Solid State Drive)

Connectivity: 1x USB 3.0 required

OS: Windows 7, 10, and 11 Pro

RECOMMENDED SYSTEM REQUIREMENTS:

CPU: Intel i7-7700 or better (8 CPU Cores or more.)

Graphics: Nvidia 4 GB Dedicated Video RAM or higher

RAM: 16GB DDR4 or higher

Storage: 1TB SSD (Solid State Drive) or larger

Connectivity: 1x USB 3.0 required

OS: Windows 7, 10, and 11 Pro

Created by: Brian Kinard
Nov 13 2021



Specialized in dental desktop scanner