


**MESTRA®**  
 POLIMERIZADORA M-9  
 POLYMERIZING M-9  
 POLYMERISEUSE M-9  
 R-080400



**MESTRA®**  
POLIMERIZADORA M-9  
**POLYMERIZING M-9**  
POLYMERISEUSE M-9  
R-080400

Nº	Ref.	DESCRIPCIÓN	DESCRIPTION
1A	080400-01	Caja de acero	Stell box
1B	080400-32	Tapa superior	Upper lid
1C	080400-33	Tapa inferior	Lower lid
1D	080400-34	Bandeja interior	Inner salver
2	000212	Cable	Cable
	000214	Cable (América)	Cable (America)
3	080010-05	Pasa cables	Cable bushing
4	000302	Cerámica	Ceramics
5	080400-06	Control (230 V)	Control panel (230 V)
	080400-06B	Control (115 V)	Control panel (115 V)
6	080400-07	Interruptor (230 V)	Switch (230 V)
	080400-10	Interruptor (115 V)	Switch (115 V)
7	080400-08	Ventilador (230 V)	Fan (230 V)
	080400-08B	Ventilador (115 V)	Fan (115 V)
8	080400-09	Casquillo	Bushing
9	080490-31	Patas	Legs
10	080251-09	Tornillos patas	Legs screws
11	080400-12	Asa	Grip
12	080460-13	Arandela de asa	Grip washer
13	080460-65	Tornillo de asa	Grip screw
14	080400-19	Grifo de desagüe	Drain cock
15	030420-22	Tuerca termopar	Thermocouple nut
18	080400-23	Termopar	Thermocouple
19	080400-24	Resistencia (230 V)	Resistor (230 V)
	080400-24B	Resistencia (115 V)	Resistor (115 V)
20	080400-25	Junta tórica resistencia y termopar	Resistor and thermocouple O-ring
21	030420-14	Arandela resistencia y termopar	Resistor and thermocouple washer
22	080400-27	Tuerca de resistencia	Resistor nut
23	080060-24	Triac	Triac
24	080400-28	Disyuntor (230 V)	Circuit breaker (230 V)
	080060-41	Disyuntor (115 V)	Circuit breaker (115 V)
25	080400-29	Contactor	Contactor