

DualTrim™ Dry Model Trimmer



Made in USA!

High Performance Rotary Dry Model Trimming Machine That Is Convertible For Wet Trimming Operation



- **Fast and Accurate Trimming**

Dry trimming eliminates the potential ecological and accuracy problems associated with wet trimming

- **Convertible for Dry or Wet Operation**

Optional package allows the conversion to wet trimming operation for maximum versatility

- **New Exhaust Port Design**

1½" OD exhaust port promotes better capture of dry and wet grinding debris resulting in less cleaning

- **New Quick-Stick™ Abrasive Disk System**

Faster trimming and fast, "wrench-less" disk changes

- **Improved Water Spray Bar Design for Wet Trimming**

Moveable spray bar promotes improved cleaning of wheel and rotates out of the way during wheel or disk changes

- **New Noise Dampening Feet**

Quieter operation and easier installation

- **Ordering Information**

<u>Item Number</u>	<u>Description</u>
--------------------	--------------------

61850	DualTrim Rotary Dry Model Trimmer, 1/3 HP, 10" Diameter, 120 V AC, 60 Hz (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61855	DualTrim Rotary Dry Model Trimmer, 1/3 HP, 10" Diameter, 220 V AC, 60 Hz (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61860	DualTrim Rotary Dry Model Trimmer, 1/2 HP, 10" Diameter, 120 V AC, 60 Hz (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61865	DualTrim Rotary Dry Model Trimmer, 1/2 HP, 10" Diameter, 220 V AC, 60 Hz (Includes 4 Quick-Stick Abrasive Disks (1 Mounted, 3 Extra), Rubber Feet and Exhaust Hose)
61875	Water Conversion Kit For DualTrim (Includes Water Spray Bar Assembly, Water Supply and Exhaust Hoses, Exhaust Hose Adapter, and Arbor Wrench Set)
78900	Custom Safety Shield Kit (Includes Clear Plastic Shield, Metal Goose-Neck, and Base that Clamps to Motor Base)
61798	Rubber Feet For Model Trimmers, Set of 4
61965	Quick-Stick Abrasive Disks, Coarse Grit, 10" Diameter (w/ 2" Center Hole), Pkg. of 6

Note: Dry model trimming requires the use of an appropriate dust collection system (not included) for removal of the dry grinding waste.