



AEU-1000-70V

General Practice

Implant Motor



- Large, super-simple, button-less touchscreen interface
- Torque control up to 50 Ncm using optional High-Efficiency 20:1 Contra Angle
- Speed Control between 15 - 2,000 RPM using 20:1 Handpiece
- Built-in, Easy-load irrigation pump with adjustable flow
- Five customizable Presets
- Compact with modern design
- Includes Multifunction Foot Control
- Worldwide Voltage Compatible
- 2-Year Warranty
- Made in USA

866.244.2954 • www.aseptico.com