

# HU-FRIEDY EMS

## BIOFILM MANAGEMENT SOLUTIONS

### THE LATEST ADVANCEMENTS IN AIRFLOW® TECHNOLOGY



#### AIRFLOW® Prophylaxis Master | FT-229HF

The AIRFLOW Prophylaxis Master is the most advanced device for delivering Guided Biofilm Therapy (GBT) to patients. Created for intensive professional use, the device features unique ergonomics, high precision, easy maintenance and compliance.



#### AIRFLOW ONE | FT-230HF

The AIRFLOW ONE shares the same meticulously designed features of the AIRFLOW Prophylaxis Master, with AIRFLOW and PERIOFLOW technology only for superior supragingival & subgingival biofilm management.

### AIRFLOW STATIONS



#### AIRFLOW STATION FEATURING WATER DELIVERY | DW-048A

The new AIRFLOW cart converts the AIRFLOW Prophylaxis Master and AIRFLOW One into an ultra-compact prophylaxis station - no water connection needed. The AIRFLOW STATION+ comes with an autonomous and independent water supply system with 360° movement for greater maneuverability, ideal for easy storage and tight spaces.



#### AIRFLOW STATION | DW-048B

The AIRFLOW STATION transforms the AIRFLOW Prophylaxis Master and AIRFLOW One into a completely mobile 360° station. Thanks to its fully-rotating wheels and minimalist design, you will benefit from greater maneuverability in the tightest of spaces.

### AIRFLOW POWDERS



#### AIRFLOW PLUS POWDER | DV-117/A

- Extra fine grains (14 µm) made from erythritol are safe for biofilm removal on all supragingival or subgingival tooth surfaces (including implants and restorations), but also effective for light to moderate stain removal
- 120g per bottle



#### PERIO DV-110/A/HS

- Fine grains (25 µm) made from glycine are gentle on tooth surface and subgingivally for efficient, safe biofilm removal
- 120g per bottle



#### CLASSIC Lemon DV-1001



#### CLASSIC Mint DV-1010



#### CLASSIC Cherry DV-1006



#### CLASSIC Neutral DV-1009

- Made from sodium bicarbonate for efficient removal of heavy stains and biofilm
- Smaller average grain size (40 µm) for increased patient comfort compared to other sodium bicarbonate powders
- 300g per bottle

### AIRFLOW HANDPIECES AND NOZZLES



#### AIRFLOW handy 3.0 | FT-220#HF/A/00x

- Easily connects to dental unit for quick setup in any operator
- Lightweight for optimal clinician comfort
- Sleek design for easy vision of treatment area
- For supragingival applications and the subgingival border area
- [x: 1 = Midwest; 2=Kavo; 3=W&H; 4=NSK; 5=Bien Air; 6/12=Sirona]



#### AIRFLOW handy 3.0 PREMIUM | FT-225#HF/A/00x

- Easily connects to dental unit for quick setup in any operator
- Lightweight for optimal clinician comfort
- Sleek design for easy vision of treatment area
- For supragingival applications and the subgingival border area
- Includes 1 PERIO handpiece and 1 PLUS handpiece
- [x: 1 = Midwest; 2=Kavo; 3=W&H; 4=NSK; 5=Bien Air; 6/12=Sirona]



#### AIRFLOW handpiece | EL-308#HF/A Handy 3.0 handpiece | EL-540#HF/A



#### PERIOFLOW handpiece | EL-354#HF/A Handy 3.0 PERIO handpiece | EL-542#HF/A



#### PERIOFLOW nozzle refills DT-476HF (coming soon)



#### Handy 3.0 Plus handpiece EL-545#HF/A

## PIEZON® TABLETOP DEVICES



**Piezon Master 700**  
| FT-194#HF/A/001

- Intelligent feedback control technology senses when more power is needed to remove deposit and adjusts accordingly
- Piezon LED handpiece with circular light emission increases vision of the treatment site



**Piezon® 250**  
FT-224#HF/A/001

- Large, ergonomic interface with fingertip control allows the clinician to easily select or change the desired power setting during use
- Illuminated handpiece remains lit for 20 seconds after the foot control is released
- Instantaneous tip amplitude adjustment offers smooth treatment and ultimate patient comfort
- Exceptionally low noise level during operation provides a more peaceful patient experience



**Piezon 150**  
FT-223#HF/A/001

## PIEZON INTEGRATED DEVICES



**Piezon 707 Built-in Kit**  
FS-435HF, FS-446HF (A-dec 541 only)

- Advanced technology of the Piezon Master 700 now available for integration into your A-dec® dental unit
- LED handpiece included

**Piezon Built-In Kit | FS-321#A/001**

- The effective Piezon technology conveniently built-in to your dental unit
- Compatible with most dental units

## PIEZON HANDPIECE



**Piezon LED handpiece**  
| EN-060#HF/A

## PIEZON INSTRUMENTS



**Instrument A**  
DS-001HF/HF

For moderate to heavy supragingival deposits



**Instrument P**  
DS-011HF/HF

For light to moderate supragingival deposits and subgingival use up to 4mm where access permits



**Instrument PS**  
DS-016HF/HF

For light to moderate subgingival deposits



To learn more visit [Hu-Friedy.com/HFEMS](http://Hu-Friedy.com/HFEMS)

©2020 Hu-Friedy Mfg. Co., LLC. All rights reserved. A-dec is a registered trademark of A-dec, Inc. AIRFLOW® is a registered trademark of Ferton Holding S.A. CORPORATION SWITZERLAND. [557\_2] HF-3000/0120

**Hu-Friedy**