



SEM picture of SphereTEC™ spherical prepolymerized fillers. Their shape and size distribution create a ball bearing-like effect that enables excellent adaptation and sculptability.



The micro-granulated structure of SphereTEC™ filler is able to bind more free resin compared to usual fillers. Therefore TPH Spectra® ST composite does not stick to hand instruments.

# TPH Spectra® ST universal composite with SphereTEC™ Technology.

Made with what matters.

Magnification of a 15 micron diameter SphereTEC™ filler particle

The key improvements in TPH Spectra® ST composite were made possible by SphereTEC™ filler technology, our proprietary method of manufacturing microsized, well-defined spherical superstructures, comprised of submicron glass. The SphereTEC™ fillers' morphology, particle size distribution and surface microstructure deliver the benefits that really matter to dentists:




## OPTIMIZED HANDLING

for easy placement and shaping



## ENHANCED CHAMELEON BLENDING ABILITY

for simplified shade selection and matching



## IMPROVED DURABILITY

with high polish and stain resistance for lasting esthetics

# TPH Spectra® ST

Universal Composite Restorative

## Ordering Information

### INTRODUCTORY KIT LV

642000

40 LV Compules® Tips (Blue) Assorted - (0.25g each)  
 5 Shade A1  
 10 shade A2  
 10 Shade A3  
 5 shade A3.5  
 5 shade A4  
 5 shade BW  
 Compules Tips Gun  
 Shade Guide  
 i-shade label



### INTRODUCTORY KIT HV

642004

40 HV Compules® Tips (Gray) Assorted - (0.25g each)  
 5 Shade A1  
 10 shade A2  
 10 Shade A3  
 5 shade A3.5  
 5 shade A4  
 5 shade BW  
 Compules Tips Gun  
 Shade Guide  
 i-shade label



### LV Refills

Shades

A1 (= VITA® shades A1; B1; C1)  
 A2 (= VITA® shades A2; B2; D2)  
 A3 (= VITA® shades A3; C2; D3; D4)  
 A3.5 (= VITA® shades A3.5; B3; B4; C3)  
 A4 (= VITA® shades A4; C4)  
 BW (Bleach White)

### Compules® Tips

(20 x 0.25g)

642220  
 642221  
 642222  
 642223  
 642224  
 642225

### Syringe

(1 x 3g)

642460  
 642461  
 642462  
 642463  
 642464  
 642465



### HV Refills

Shades

A1 (= VITA® shades A1; B1; C1)  
 A2 (= VITA® shades A2; B2; D2)  
 A3 (= VITA® shades A3; C2; D3; D4)  
 A3.5 (= VITA® shades A3.5; B3; B4; C3)  
 A4 (= VITA® shades A4; C4)  
 BW (Bleach White)

### Compules® Tips

(20 x 0.25g)

642231  
 642232  
 642233  
 642234  
 642235  
 642236

### Syringe

(1 x 3g)

642470  
 642471  
 642472  
 642473  
 642474  
 642475



### Shade Guide Refill

642013

1 Shade Guide Holder  
 6 shade guide tabs (A1, A2, A3, A3.5, A4, BW)



### COMPULES TIPS GUN

635105

1 Compules Tips Gun



NEW  
with SphereTEC™  
Filler Technology

THE NEWEST BREAKTHROUGH  
IN COMPOSITE TECHNOLOGY  
BEGAN WITH A QUESTION:

**WHAT MATTERS TO YOU?**

Magnification of a 15 micron diameter SphereTEC™ filler particle

TPH Spectra® ST  
Universal Composite Restorative

Dentsply  
Sirona

A dental handpiece is shown in the upper right, angled downwards towards a tooth model on the left. The handpiece is dark and metallic. The tooth model is a realistic, light-colored dental specimen. The background is a light, neutral color. Several semi-transparent blue circles of varying sizes are overlaid on the tooth model and the handpiece. In the bottom right corner, there are several semi-transparent white circles of varying sizes, representing the filler particles mentioned in the text.

We didn't just add  
SphereTEC™ filler  
to TPH Spectra®  
composite, we  
added more of  
**what matters**  
**to you.**

## TPH Spectra® ST

Universal Composite Restorative

When we set out to improve TPH Spectra® universal composite, we asked clinicians what matters most to them. Certainly, it needs to be durable. But could the material handle better to save you time and frustration during placement? And what if shade selection was simplified, minimizing stock and creating efficiencies during the procedure? Could this composite offer outstanding stain resistance to keep patients satisfied?

**Introducing TPH Spectra® ST, composite made to deliver what matters to you.**

# efficiency matters

74% of dentists complain about handling properties of composites such as adaptability and stickiness.<sup>2</sup>

## Versatile handling, made even easier to work with.

Thanks to new SphereTEC™ filler technology, TPH Spectra® ST composite offers optimized handling characteristics in both a high and low viscosity to accommodate your preference and technique for placement efficiency.

- Doesn't stick to hand instruments
- Adapts easily to cavity surfaces
- Easy to sculpt and shape
- Excellent slump resistance



# esthetics matter

Dentists spend up to 28% of Class II chair time on shade selection and finishing and polishing.<sup>1</sup>

## Exceptional esthetics, made even easier to create.

SphereTEC™ filler brings an enhanced chameleon blending ability to TPH Spectra® ST composite for simplified shade matching. It also enables faster finishing and polishes to an excellent luster.

- Five CLOUD shades cover the full VITA® range<sup>3</sup>
- Save money and shelf space by reducing inventory
- Polish to a high gloss, faster with the Enhance® Finishing and Enhance® PoGo® Polishing systems



## FIVE CLOUD SHADES

VITA® Shades<sup>1</sup>

A1			A2			A3				A3.5				A4		BW (BLEACH WHITE)
A1	B1	C1	A2	B2	D2	A3	C2	D3	D4	A3.5	B3	B4	C3	A4	C4	
A1	B1	C1	A2	B2	D2	A3	C2	D3	D4	A3.5	B3	B4	C3	A4	C4	

# longevity matters

## Esthetic results made even more durable.

With 24 years and over 174 million restorations placed, the TPH® brand has a proven record of long-term clinical success. New SphereTEC™ filler takes that performance to the next level.

- Excellent stain resistance
- 60% improved wear resistance

57% of dentists frequently see composite staining at recall appointments.<sup>5</sup>

80% of patients are reportedly aware of color differences between restored and adjacent natural teeth.<sup>4</sup>



After 24-Hour Staining Test<sup>1</sup>

TPH Spectra® ST Composite with its new SphereTEC™ filler technology is made for what matters to you.