

6700M Electro-mechanic Chair

Benefits of the 6700M Electro-Mechanic Operation

- Precision: Precise and smooth movements with no hard start/stops
- Easy Diagnostic: Simple feedback loops make it easier to monitor and measure the health of the system
- Energy Efficient: Less power usage than a hydraulic pump which means reduced maintenance cost and reduced operating cost
- Simple: Less complex direct drive
- Clean: No leaks and no liquids involved
- Quiet: Almost soundless operation
- Safe: No chance of high-pressure oil leak
- Durable: Millions of cycles possible
- Weight Capacity: Equal to hydraulic pumps




Features and Functions

	6700M Chair
Construction	Steel construction
Auto Exit and Pre-Programmable Positions	Auto exit four pre-programmable positions Last position
Operations/Movement	Electro-mechanic base & back
Controls	Moveable foot control and side chair touch pads
Swivel Base	Not applicable
Armrests	Swivel out arms
Headrest	Double articulating knob style headrest

Internal Components

	All Electro-mechanic Models
Access	All components are readily accessible
Operation	Quiet German motors
PC Board	Simple diagnostics with remote diagnostics and USB access, easily programmable features



800.275.3368 | Tel: 954.730.3636 | Fax 954.730.3602
1280 SW 27th Ave | Pompano Beach, FL 33069
www.summitdental.com

