

Cyclone Dust Collector

Designed for trouble-free performance and heavy duty operation. Quiet, portable, and compact, the dust collector is ideal for under the counter operation.

Features

- Very quiet 59 dB noise level.
- One HP motor, 112 cubic feet per minute (CFM) suction power
- Economical, disposable paper primary filter, easily and quickly changed. 2-ply filter retains 1 micron particles.
- Permanent cloth secondary filter.
- Supplied with one 6' long x 2-1/2" ID (1.8 m x 63.5 mm) vacuum hose and four (4) extra disposable paper dust bags.
- Optional remote foot switch (item X131) CDC1 model, included with CDC2 model.
- Uses 1-1/2", 2-1/2", or 3" I.D. hose for adaptability to nearly every dust collection devise.



Specifications	Item CDC1
Unit size: 12-1/2" x 12-1/2" x 17" high (32 cm x 32 cm x 43 cm)	
115 volt, 60 hz standard	
230 volt, 50/60 hz optional for export	
35 lbs. (16 kg) net weight,	
42 lbs. (19 kg) shipping weight	

The Ray Foster Cyclone Dust Collector is the perfect compliment to your lab environment and can be directly connected to the following Ray Foster equipment for automatic on/off operation of the Cyclone Dust Collector.

AG03 and AG05 Alloy Grinders

MTDI Dry Model Trimmer

MR20 Model Arch Trimmer



The Deluxe Cyclone Dust Collector

featuers a two stage by-pass motor, suitable for use with abrasive blaster, and includes a remote foot switch and mobile cart.

Specifications	Item CDC2
Unit size: 12-1/2" x 12-1/2" x 20" high (32 cm x 32 cm x 51 cm)	
115 volt, 60 hz standard	
230 volt, 50/60 hz optional for export	
42 lbs. (19 kg) net weight,	
50 lbs. (23 kg) shipping weight	



Item X210LS Splash Hood

Steel construction with stainless steel tray insert, light, shield, and flange for suction connection.
6"w. x 10"d. x 12"h. (12.5 cm x 25 cm x 30.5 cm)



Item XI48 Counter Air Scoop

Cast aluminum with screened port. Mouth measures 7" wide by 5-1/2" high (18 cm x 14 cm).



Item XI52 Mobile Cart

Four swivel casters mounted on a steel frame, will make your cyclone Dust Collector truly portable.



Item XI31 Remote Foot Switch

Heavy duty foot switch for convenient hands-free on/off control.