

Ray Foster Electric Wax Heaters



Catalog WH41

Deluxe Wax Heater with Thermostat Control

Full range of temperature adjustment from 100°F (38°C) through 250°F (120°C).

Ray Foster Electric Wax Heaters feature three compartments with two measuring 2-1/8" (54mm) X 1-3/4" (45mm) X 3/4" (19mm) deep, and one 2-1/4" (57mm) X 3-3/4" (95mm) X 3/4" (19mm) deep. They have a hinged lid to keep out contaminants, indicator lights, and an ON/OFF switch. Rugged construction and reliable components promise years of trouble free service. Wax Heaters are available in two models.



Catalog WH40

Standard Wax Heater

Single preset temperature of 160° (71°C).

www.fosterdental.com