

# Ray Foster Dental Vibrators



**Catalog DV34**  
*Variable Speed Dental Vibrator*

With the variable speed vibration control feature this vibrator is excellent for general purpose use. It is solidly constructed with a cast metal housing, 4-1/4" (10.8 cm) diameter removable plastic top, and shock absorbing rubber feet. Net weight is 5 lb, 2 oz (2.33 kg).



**Catalog DV36**  
*Heavy Duty Dental Vibrator*

Our DV36 Vibrator provides high vibration intensity with a three-speed switch. It has a removable rubber top measuring 6-1/2" (16.5 cm) diameter. Net weight is 9 lbs (4.1 kg).

**R**ay Foster offers a choice of two dental vibrators perfect for any laboratory. From the gentlest vibration to a powerful shake, these vibrators offer quiet, creep-free performance for almost any application.

[www.fosterdental.com](http://www.fosterdental.com)