

# Renfert

making work easy

## NEW POWER steamer 1

Steam cleaning unit with manual filling

The new, powerful POWER steamer 1 sets new standards in reliability and durability and impresses with its outstanding cleaning performance. Any kind of typical soiling can be easily and effectively removed. The location can be selected flexibly thanks to manual filling.

### Advantages

- Particularly durable heating system – the heating element is cast inside the floor of the boiler to protect it against water, limescale, and the chemicals used during cleaning.
- High functional reliability and dependability thanks to a real-time calcification indicator with 2 increments.
- Simple cleaning of the pressure tank thanks to an extra-large service opening.
- Smart operation plus additional functions via the Wi-Fi interface and Renfert CONNECT stick (optional accessory)\*.

\*Renfert CONNECT is not available in all regions.

### Details

- Continuous and effective cleaning performance thanks to the 2,000 W and 4.5 bar working pressure in combination with innovative steam nozzle technology.
- Flexible location selection thanks to manual filling.
- Safe and convenient steam cleaning thanks to an ergonomic handpiece (the housing safely dissipates electrostatic charge).
- Easy filling of water and descaling solution thanks to handy, hopper-shaped fillers.
- Wall mounting possible using an optional wall bracket.



tool  
included



easy use



compact



smart  
control



Innovative steam nozzle technology for effective cleaning performance

### Technical data

Permissible mains voltage	120 V
Permissible mains frequency	50/60 Hz
Pressure vessel volume	1.06 gal
Volume (Recommended fill quantity)	0.74 gal
Volume (Maximum fill quantity)	0.79 gal
Working pressure	65.3 psi
Dimensions (W x H x D)	14.6 x 16.6 x 13.8"
Weight (empty)	~20.9 lbs

### Ordering information

POWER steamer 1 120V	No. 18451000
Wall bracketed	No. 18450200
POWER steamer descaler	No. 18450100

# Renfert

making work easy

## NEW POWER steamer 2

Steam cleaning unit with automatic filling

The new, powerful POWER steamer 2 sets new standards in reliability and durability and impresses with its outstanding cleaning performance. Any kind of typical soiling can be easily and effectively removed. Using the mains water connection and integrated pump, the pressure tank is filled automatically and ensures availability of steam at all times.

### Advantages

- Particularly durable heating system – the heating element is cast inside the floor of the boiler to protect it against water, limescale, and the chemicals used during cleaning.
- High functional reliability and dependability thanks to a real-time calcification indicator with 2 increments.
- Simple cleaning of the pressure tank thanks to an extra-large service opening.
- Smart operation plus additional functions via the Wi-Fi interface and Renfert CONNECT stick (optional accessory)\*.

\*Renfert CONNECT is not available in all regions.

### Details

- Continuous and effective cleaning performance thanks to the 2,000 W and 4.5 bar working pressure in combination with innovative steam nozzle technology.
- Availability of steam at all times – even during heavy-duty continuous use – thanks to automatic filling using the mains water connection and integrated pump.
- Optimal control through display of the current working pressure using a manometer.



tool included



easy use



compact



smart control



Practical storage of cleaning items

### Technical data

Permissible mains voltage	120 V
Permissible mains frequency	50/60 Hz
Pressure vessel volume	1.06 gal
Volume (Recommended fill quantity)	0.74 gal
Volume (Maximum fill quantity)	0.79 gal
Working pressure	65.3 psi
Connection pressure	14.5–72.5 psi
Dimensions (W x H x D)	14.6 x 16.6 x 13.8"
Weight	~22.3 lbs

### Ordering information

POWER steamer 2 120V	No. 18461000
Wall bracketed	No. 18450200
POWER steamer descaler	No. 18450100
POWER steamer water softener	No. 18460100