

# P42 Grinding machine

Intelligent ,Easy and Precise!



01 seeking for better aesthetic results

Support ultra-thin veneers (0.15mm) processing

02 save time and cost

Support glass ceramic grinding with clear water

03 A-axis can swing  $\pm 20^\circ$

make undercut process better

04 Noise reduction design

Perfect choice for clinics & in-house labs

## P42 Grinding machine



Model	P42
Dimension	L*W*H: 560x400x565 (mm)
Weight	80KG
Input Voltage	AC 220V±10%、50~60Hz
Product Power	0.8KW
Spindle Power	0.55 KW
Axis	4
Milling Range	X/Y/Z: 110/75/50MM A: ±360°
Mode	Wet Milling
Max Spindle Speed	60,000 rpm
Number of Tools	6
Certifications	CE、Fcc、Rohs
Compatible CAM Software	UPCAM、Millbox
Milling Materials	Glass Ceramic、Titanium
Operation Temperature	10 °C - 50 °C
Operation Humidity	< 93% RH (no condensation)

**UP3D | CONTACT US**

Web: [www.up3ds.com](http://www.up3ds.com)

Email: [info@up3d.cn](mailto:info@up3d.cn)

Tel: +86-755-26983202

