



ITEM NO. 72670 12 PACK

LOT # m029010

.040" (1mm) RED DENTAL EVA FOR
BUILDUP OF MOUTHGUARDS

ROUND

-FOR ATHLETIC AND RECREATIONAL USE-

(Completed Athletic Mouthguards should be at least .160" in thickness)
(125mm Diameter Disk)

Heat to allow slight sag at center. Do not overheat.
Store in a cool area.

NO PVC

NO BPA (used in UV cured resins)
NO LATEX / NO POLYCARBONATE
NO CHINA CONTENT

Made in USA

Manufactured by: Buffalo Dental Mfg., Syosset, NY 11791

BUFFALO

424958

PolyShoK .155" (4mm) Black Sports M

72670

www.UsDentalDepot.com